


Dr. Maisch

Any Column, Any Size, Any Media



REPROSHELL®

Core-Shell columns

MADE BY DR. MAISCH



CONTENT

- P 4 - REPROSHELL - CORE-SHELL-PARTICLES
- P 5 - REPROSHELL - VAN DEEMTER EQUATION
- REPROSHELL - VAN DEEMTER PLOT
- P 6 - A-TERM (EDDY DIFFUSION, MULTI-PATH EFFECT)
- B-TERM (LONGITUDINAL DIFFUSION)
- C-TERM (MASS TRANSFER)
- P 7 - COMPARISON CORE-SHELL vs. FULLY POROUS PARTICLES
- P 8 - REPROSHELL - AVAILABLE PHASES
- P 9 - 11 - REPROSHELL - EFFICIENCY
- P 12 - 15 - AVAILABLE COLUMNS



**REPROSHELL
MADE BY DR. MAISCH**


From one of the largest **H**igh-**P**erformance
Liquid **C**hromatography (HPLC) - Column
Manufacturers in Europe.

REPROSHELL

Core - Shell particles




Core-Shell particles consist of a solid (non-porous) core surrounded by a porous layer (shell). They exhibit the same back pressure but higher efficiency compared to fully porous particles of the same particle size.




REPROSHELL

ReproShell – van Deemter Equation



van Deemter Plot




REPROSHELL

A-TERM (EDDY DIFFUSION, MULTI-PATH EFFECT)




Core-Shell particles have a tighter particle size distribution which leads to a more uniform, organised bed structure. This reduces the impact of the multi-path effect.

B-Term (Longitudinal Diffusion)



Fully porous particles:


Diffusion takes place within the porous particle as well as in the interstitial space



Core-Shell particles:

Non-porous core blocks longitudinal diffusion


C-Term (Mass transfer)



The non-porous core leads to a shorter diffusion path

REPROSHELL

COMPARISON CORE-SHELL VS. FULLY POROUS PARTICLES




REPROSHELL

ReproShell – Available phases

Particle sizes 2.7 μm and 5 μm



UPLC PERFORMANCE WITH STANDARD HPLC INSTRUMENT
AND BACKPRESSURE
VARIOUS SELECTIVITIES

- ODS-1 (highly retentive C18)
- ODS-3 (trifunctional C18)
- C8
- Biphenyl
- Phenylhexyl
- PFP (Penta fluoro phenyl)
- Silica
- HILIC




REPROSHELL - EFFICIENCY



ReproShell ODS - 1
2.7 μm





ReproShell ODS - 1
5 μm



ReproShell ODS - 3
2.7 μm





ReproShell ODS - 3
3 μm





REPROSHELL - EFFICIENCY



ReproShell Biphenyl,
2.7 μm





ReproShell Biphenyl,
5 μm



ReproShell C8,
2.7 μm




ReproShell C8,
5 μm



REPROSHELL - EFFICIENCY

HPLC-UV quantitation of 11 cannabinoids



- | | | |
|------------|----------|----------|
| 1 - THCA | 5 - CBGA | 9 - CBDV |
| 2 - D9 THC | 6 - CBG | 10 - CBN |
| 3 - D8 THC | 7 - CBDA | 11 - CBC |
| 4 - THCV | 8 - CBD | |

Column:	Reproshell ODS-1, 2.7 μm um, 150 x 4,6 mm
Mobile Phase A:	Acetonitrile w / 0.1 % formic acid
Mobile Phase B:	Water w / 0.1 % formic acid
Mobile Phase Composition:	75/25 - A/B - Isocratic
Flow rate:	2 mL/min
Back pressure:	220-230 bar
Detection:	228 nm
Run Time:	7 minutes
Reference Standard:	Cayman Chemical, Prod No. 32842

Chromatogramme and method courtesy of www.seedlessanalytical.com

AVAILABLE COLUMNS

Coreshell RP-Phase with nonporous core and porous shell

Particle size 5µm	20 x 2 mm	30 x 2 mm	50 x 2 mm
ReproShell ODS-1	cs15.91.s0202	cs15.91.s0302	cs15.91.s0502
ReproShell ODS-3	cs15.93.s0202	cs15.93.s0302	cs15.93.s0502
ReproShell C8	cs15.8e.s0202	cs15.8e.s0302	cs15.8e.s0502
ReproShell Biphenyl	cs15.bpe.s0202	cs15.bpe.s0302	cs15.bpe.s0502
ReproShell Phenylhexyl	cs15.ph.s0202	cs15.ph.s0302	cs15.ph.s0502
ReproShell PFP	cs15.pfp.s0202	cs15.pfp.s0302	cs15.pfp.s0502
ReproShell Silica	cs15.00.s0202	cs15.00.s0302	cs15.00.s0502

Particle size 2,7 µm	20 x 3 mm	30 x 3 mm	50 x 3 mm
ReproShell ODS-1	cs27.91.s0203	cs27.91.s0303	cs27.91.s0503
ReproShell ODS-3	cs27.93.s0203	cs27.93.s0303	cs27.93.s0503
ReproShell C8	cs27.8e.s0203	cs27.8e.s0303	cs27.8e.s0503
ReproShell Biphenyl	cs27.bpe.s0203	cs27.bpe.s0303	cs27.bpe.s0503
ReproShell Phenylhexyl	cs27.ph.s0203	cs27.ph.s0303	cs27.ph.s0503
ReproShell PFP	cs27.pfp.s0203	cs27.pfp.s0303	cs27.pfp.s0503
ReproShell Silica	cs27.00.s0203	cs27.00.s0303	cs27.00.s0503

AVAILABLE COLUMNS

75 x 2 mm

cs15.91.s0702
cs15.93.s0702
cs15.8e.s0702
cs15.bpe.s0702
cs15.ph.s0702
cs15.pfp.s0702
cs15.00.s0702

100 x 2 mm

cs15.91.s1002
cs15.93.s1002
cs15.8e.s1002
cs15.bpe.s1002
cs15.ph.s1002
cs15.pfp.s1002
cs15.00.s1002

125 x 2 mm

cs15.91.s1202
cs15.93.s1202
cs15.8e.s1202
cs15.bpe.s1202
cs15.ph.s1202
cs15.pfp.s1202
cs15.00.s1202

150 x 2 mm

cs15.91.s1502
cs15.93.s1502
cs15.8e.s1502
cs15.bpe.s1502
cs15.ph.s1502
cs15.pfp.s1502
cs15.00.s1502

5 x 2 mm 3/pck

cs15.91.v0002
cs15.93.v0002
cs15.8e.v0002
cs15.bpe.v0002
cs15.ph.v0002
cs15.pfp.v0002
cs15.00.v0002

75 x 3 mm

cs27.91.s0703
cs27.93.s0703
cs27.8e.s0703
cs27.bpe.s0703
cs27.ph.s0703
cs27.pfp.s0703
cs27.00.s0703

100 x 3 mm

cs27.91.s1003
cs27.93.s1003
cs27.8e.s1003
cs27.bpe.s1003
cs27.ph.s1003
cs27.pfp.s1003
cs27.00.s1003

125 x 3 mm

cs27.91.s1203
cs27.93.s1203
cs27.8e.s1203
cs27.bpe.s1203
cs27.ph.s1203
cs27.pfp.s1203
cs27.00.s1203

150 x 3 mm

cs27.91.s1503
cs27.93.s1503
cs27.8e.s1503
cs27.bpe.s1503
cs27.ph.s1503
cs27.pfp.s1503
cs27.00.s1503

5 x 3 mm 3/pck

cs27.91.v0003
cs27.93.v0003
cs27.8e.v0003
cs27.bpe.v0003
cs27.ph.v0003
cs27.pfp.v0003
cs27.00.v0003

AVAILABLE COLUMNS

Particle size 2,7 µm	20 x 4 mm	30 x 4 mm	50 x 4 mm
ReproShell ODS-1	cs27.91.s0204	cs27.91.s0304	cs27.91.s0504
ReproShell ODS-3	cs27.93.s0204	cs27.93.s0304	cs27.93.s0504
ReproShell C8	cs27.8e.s0204	cs27.8e.s0304	cs27.8e.s0504
ReproShell Biphenyl	cs27.bpe.s0204	cs27.bpe.s0304	cs27.bpe.s0504
ReproShell Phenylhexyl	cs27.ph.s0204	cs27.ph.s0304	cs27.ph.s0504
ReproShell PFP	cs27.pfp.s0204	cs27.pfp.s0304	cs27.pfp.s0504
ReproShell Silica	cs27.00.s0204	cs27.00.s0304	cs27.00.s0504

Particle size 2,7 µm	20 x 4.6 mm	30 x 4.6 mm	50 x 4.6 mm
ReproShell ODS-1	cs27.91.s0246	cs27.91.s0346	cs27.91.s0546
ReproShell ODS-3	cs27.93.s0246	cs27.93.s0346	cs27.93.s0546
ReproShell C8	cs27.8e.s0246	cs27.8e.s0346	cs27.8e.s0546
ReproShell Biphenyl	cs27.bpe.s0246	cs27.bpe.s0346	cs27.bpe.s0546
ReproShell Phenylhexyl	cs27.ph.s0246	cs27.ph.s0346	cs27.ph.s0546
ReproShell PFP	cs27.pfp.s0246	cs27.pfp.s0346	cs27.pfp.s0546
ReproShell Silica	cs27.00.s0246	cs27.00.s0346	cs27.00.s0546

AVAILABLE COLUMNS

75 x 4 mm

cs27.91.s0704
cs27.93.s0704
cs27.8e.s0704
cs27.bpe.s0704
cs27.ph.s0704
cs27.pfp.s0704
cs27.00.s0704

100 x 4 mm

cs27.91.s1004
cs27.93.s1004
cs27.8e.s1004
cs27.bpe.s1004
cs27.ph.s1004
cs27.pfp.s1004
cs27.00.s1004

125 x 4 mm

cs27.91.s1204
cs27.93.s1204
cs27.8e.s1204
cs27.bpe.s1204
cs27.ph.s1204
cs27.pfp.s1204
cs27.00.s1204

150 x 4 mm

cs27.91.s1504
cs27.93.s1504
cs27.8e.s1504
cs27.bpe.s1504
cs27.ph.s1504
cs27.pfp.s1504
cs27.00.s1504

5 x 4 mm 3/pck

cs27.91.v0004
cs27.93.v0004
cs27.8e.v0004
cs27.bpe.v0004
cs27.ph.v0004
cs27.pfp.v0004
cs27.00.v0004

75 x 4.6 mm

cs27.91.s0746
cs27.93.s0746
cs27.8e.s0746
cs27.bpe.s0746
cs27.ph.s0746
cs27.pfp.s0746
cs27.00.s0746

100 x 4.6 mm

cs27.91.s1046
cs27.93.s1046
cs27.8e.s1046
cs27.bpe.s1046
cs27.ph.s1046
cs27.pfp.s1046
cs27.00.s1046

125 x 4.6 mm

cs27.91.s1246
cs27.93.s1246
cs27.8e.s1246
cs27.bpe.s1246
cs27.ph.s1246
cs27.pfp.s1246
cs27.00.s1246

150 x 4.6 mm

cs27.91.s1546
cs27.93.s1546
cs27.8e.s1546
cs27.bpe.s1546
cs27.ph.s1546
cs27.pfp.s1546
cs27.00.s1546

5 x 4.6 mm 3/pck

cs27.91.v0046
cs27.93.v0046
cs27.8e.v0046
cs27.bpe.v0046
cs27.ph.v0046
cs27.pfp.v0046
cs27.00.v0046

Dr. Maisch

Any Column, Any Size, Any Media

Distributor:



Dr. Maisch HPLC GmbH
Beim Brückle 14
D-72119 Ammerbuch
T: +49 7073 50357
F: +49 7073 4216
www.Dr-Maisch.com
www.MODcol.com
E-Mail: info@dr-maisch.com



PDF brochure
for download