

Organizers

HPLC Columns and Accessories Organizers

- Holds 10 HPLC columns
- 17 compartments
- PVC construction
- 12" H x 24" W x 6" D
(30.48cm H x 60.96cm W x 15.24cm D)



HPLC Columns Organizers

- 30 HPLC columns
- Wall-mounts with two screws
- PVC construction with acrylic doors
- 17-3/8" H x 15" W x 2-7/8" D
(44.18cm H x 38.1cm W x 7.3cm D)



Plastic Box with Removable Dividers

- 13 removable dividers
- Ideal for fittings and bore-cut tubing



Lab Organizer with Adjustable Bins

- 14 adjustable dividers, 4 shelves
- PVC construction
- 20" H x 24" W x 10.5" D
(50.8cm H x 60.9cm W x 26.7cm D)



GC Column and Accessory Organizers

- 5-Column Organizer:
13" H x 10.5" W x 14" L
(33cm H x 26.5cm W x 35.5cm L)
- 10-Column Organizer:
12.75" H x 10.5" W x 25" L
(32cm H x 26.5cm W x 63.5cm L)



Organizer Box

- 9-compartment organizer includes eight 1-3/8" x 1-5/8" (3.5cm x 4.13cm) slots and one 1-1/8" x 3-1/4" (2.8cm x 8.2cm) slot
- 18-compartment organizer includes twelve 2" x 2" (5cm x 5cm) slots and six 2" x 4-3/16" (5cm x 10.6cm) slots



Organizers

Description	Part No.
HPLC Supplies Organizer	5238
HPLC Column Storage Unit	5240
Plastic Storage Box with Removable Dividers	803880
Lab Organizer with Adjustable Bins	5242
GC Column and Accessory Organizers	
GC Column Organizer, 5-Column	6390
GC Column Organizer, 10-Column	6392
Organizer Box	
9-Compartment Box	6264
18-Compartment Box	6266

related products


Looking for vial organizers and storage racks?
See pages 380 and 381.



Tools


Magnifiers

- For inspecting capillary column ends



Deluxe Torch Kit

- Ideal for sealing the ends of capillary columns
- Kit contains the torch, six Micronox cylinders, three butane cylinders, six silver soldering rods, flux tube, two flame tips, one butane tip, one spark lighter, and complete instructions



Capillary Tool Kit

- Kit includes diamond scribe, two open end wrenches, PTFE-coated forceps, ruler, ferrule pick, 10' high temperature string, 20X magnifier, and pipe cleaner



23-Piece Tool Kit in Zippered Case

- Kit includes seven screwdrivers, adjustable wrench, two pliers, wire stripper, knife, alignment tool, stainless rule, hex key set, scissors, two flexible files, burnisher, miniature soldering iron, solder aid, coil of solder, and desoldering braid



Adjustable Wrench

- Chrome-Plated



Open-End Wrench

- Chrome-Plated



Handy-Lok Fitting Wrench

- Ideal for connections on injection valves



9-Way Wrench

- For use with gas cylinders



Miniature Files, Fine Cut

- Set includes six different-shaped, 5-1/2" long files
- Made of hard chrome steel suitable for smoothing and deburring



Miniature Reamers

- Set includes six different sizes: 0.028"/0.018", 0.040"/0.026", 0.048"/0.032", 0.080"/0.050", 0.085"/0.053", 0.104"/0.062"



Flame Checker

- Pivoting mirror head fogs when flame is lit



Tools

Description	Qty.	Part No.
Magnifier, 20X	ea	3163
Magnifier, Lighted (includes battery)	ea	3196
Deluxe Torch Kit		14606
Micronox Cylinders	6/pk	14602
Butane Cylinders	6/pk	14603
Piercing Knob, Butane	ea	14604
Piercing Knob, Micronox	ea	14605
Capillary Tool Kit	ea	16879
23-Piece Tool Kit and Case	ea	16533
Adjustable Wrench		
4" Long—1/2" Capacity	ea	16549
6" Long—3/4" Capacity	ea	3211
8" Long—7/8" Capacity	ea	16522
Open-End Wrench		
1/4" x 5/16"—Fits 1/16" Fitting Nuts	ea	1998
7/16" x 1/2"—Fits 1/8" and 3/16" Fitting Nuts	ea	1999
9/16" x 5/8"—Fits 1/4" and 3/8" Fitting Nuts	ea	2000
Handy-Lok Fitting Wrench	ea	35751
9-Way Cylinder Wrench	ea	412933
Miniature Files, Fine Cut	6/pk	7237
Miniature Reamers	6/pk	7330
Flame Checker	ea	7047