

HPLC Starter Kit 3 V2.0

SafetyCap 2.0

GL 45 Thread

Universal!



**THE System
for Solvent Safety**

**Initial equipment for your HPLC system -
Safe extraction of solvents.**



- Price advantage compared to ordering piece by piece.
- Suitable for all current HPLC systems.
- Unused connectors can be closed by the blind screws delivered with the set. This makes the system universally applicable, even when less connectors are needed for a time.



Part No.	Description	Thread size	Connections for Ø 1.6 / 2.3 / 3.2 mm OD tubing	Service life (in months)	Colour	Sales Unit ----- incl.
399 202	HPLC Starter Kit 3, V2.0					1
incl.	SafetyCap III, V2.0	GL 45	3		blue/white	4x
incl.	Air valve for SafetyCaps, V2.0	UNF 1/4" 28G		6 months each	red	4x
incl.	Fitting 1.6 mm	UNF 1/4" 28G			green	12x
incl.	Fitting 2.3 mm	UNF 1/4" 28G			violet	12x
incl.	Fitting 3.2 mm	UNF 1/4" 28G			blue	12x
incl.	Blind plug	UNF 1/4" 28G			translucent	8x
incl.	Insertion Label / Title Block				white	8x

www.scat-europe.com/new-generation-2-0

Phone: +49 (0) 61 05 / 30 55 86 - 0 / e-Mail: info@scat-europe.com / Internet: www.scat-europe.com

S.C.A.T. Europe GmbH / Waldecker Straße 7 / 64546 Mörfelden-Walldorf / Germany

#200144 / 15072020