

**Agilent J&W LTM II Low Thermal Mass Column Modules  
for 7890A/B Series GC Systems**

Phase	ID (mm)	Length (m)	Film (µm)	Part No.
HP-5ms	0.18	20	0.18	19091S-577LTM
		12	0.33	19091S-101LTM
		25	0.33	19091S-102LTM
	0.25	15	0.10	19091S-331LTM
			0.25	19091S-431LTM
		30	0.25	19091S-433LTM
	0.32	10	0.50	19091S-111LTM
		30	0.25	19091S-413LTM
HP-5ms Ultra Inert	0.18	20	0.18	19091S-577UILTM
		15	0.25	19091S-431UILTM
			0.25	19091S-433UILTM
			0.50	19091S-133UILTM
	0.32	30	0.25	19091S-413UILTM
			1.00	19091S-213UILTM
		30	1.00	19091S-213UILTM
HP-88	0.25	30	0.20	112-8837LTM
HP-Fast Residual Solvent	0.53	30	1.00	19095V-420LTM
HP-FFAP	0.20	25	0.33	19091F-102LTM
		30	0.25	19091F-433LTM
	0.32	25	0.50	19091F-112LTM
		10	1.00	19095F-121LTM
			30	1.00
HP-INNOWax	0.18	20	0.18	19091N-577LTM
		25	0.20	19091N-102LTM
	0.25	5	0.15	19091N-030LTM
		30	0.25	19091N-133LTM
	0.32	30	0.15	19091N-013LTM
		30	1.00	19095N-123LTM

(Continued)

