













## ZORBAX 80 Å StableBond

Hardware Beschreibung	LC Säulen- abm. (mm)	Partikel- größe (µm)	SB-C18 USP L1	SB-C8 USP L7	SB-CN USP L10	SB-C3 USP L56	SB-Phenyl USP L11	SB-Aq
<b>Standardsäulen (keine besondere Hardware erforderlich)</b>								
Narrowbore RRHT, 600 bar	2,1 x 150	1,8	820700-902	820700-906	820700-905		820700-912	
Narrowbore RRHT, 600 bar	2,1 x 100	1,8	828700-902	828700-906	828700-905		828700-912	828700-914
Narrowbore RRHT, 600 bar	2,1 x 50	1,8	827700-902	827700-906	827700-905		827700-912	827700-914
Narrowbore RRHT, 600 bar	2,1 x 30	1,8	824700-902	824700-906	824700-905		824700-912	824700-914
Narrowbore RRHT, 600 bar	2,1 x 20	1,8	826700-902	826700-906				
 UHPLC-Vorsäule, 1200 bar, 3 St./VE	2,1 x 5	1,8	821725-902	821725-904				
MicroBore RR	1,0 x 150	3,5	863600-902	863600-906	863600-905			
MicroBore RR	1,0 x 50	3,5	865600-902	865600-906				
MicroBore RR	1,0 x 30	3,5	861600-902	861600-906				
MicroBore Vorsäule, 3 St.	1,0 x 17	5	5185-5920	5185-5920				
 Vorsäulenkartuschen, 2 St.	9,4 x 15	7	820675-115	820675-115	820675-124	820675-124	820675-115	
 Vorsäulenkartuschen, 4 St.	4,6 x 12,5	5	820950-920	820950-915	820950-916	820950-922	820950-917	820950-933
 Vorsäulenkartuschen, 4 St.	2,1 x 12,5	5	821125-915	821125-915	821125-924	821125-924	821125-915	821125-933
 Vorsäulen Hardware-Kit	9,4 x 15	0	840140-901	840140-901	840140-901	840140-901	840140-901	
 Vorsäulen Hardware-Kit			820999-901	820999-901	820999-901	820999-901	820999-901	820999-901
<b>PrepHT-Kartuschensäulen (Endfitting-Kit 820400-901 erforderlich)</b>								
 PrepHT Kartusche	21,2 x 250	7	877250-102	877250-106	877250-105		877250-112	877250-114
 PrepHT Kartusche	21,2 x 150	7	877150-102	877150-106				877150-114
 PrepHT Kartusche	21,2 x 150	5	870150-902	870150-906				870150-914
 PrepHT Kartusche	21,2 x 100	5	870100-902	870100-906				870100-914
 PrepHT Kartusche	21,2 x 50	5	870050-902	870050-906				870050-914
 PrepHT Vorsäulenkartusche, 2 St.	17,0 x 7,5	5	820212-920	820212-915	820212-915		820212-915	820212-933
Vorsäulenkartuschen- Hardware			820444-901	820444-901	820444-901	820444-901	820444-901	820444-901
PrepHT Endfittings, 2 St.			820400-901	820400-901	820400-901	820400-901	820400-901	820400-901

Falls nicht anders angegeben beträgt das Drucklimit für die Säulen 400 bar.