

## Rapid Resolution HT-Säulen zur Verwendung bei hohem Druck (Maximaler Druck: 600 bar, 9000 psi)

| Hardware                                                                            | Beschreibung                             | LC Säulen-<br>abm. (mm) | Partikel-<br>größe<br>(µm) | Eclipse<br>Plus C18<br>USP L1 | Eclipse<br>Plus C8<br>USP L7 | Eclipse Plus<br>Phenyl Hexyl<br>USP L11 | Eclipse<br>PAH<br>USP L1 | Eclipse<br>XDB-C18<br>USP L1 | Eclipse<br>XDB-C8<br>USP L7 | Extend-C18<br>USP L1 |
|-------------------------------------------------------------------------------------|------------------------------------------|-------------------------|----------------------------|-------------------------------|------------------------------|-----------------------------------------|--------------------------|------------------------------|-----------------------------|----------------------|
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 150               | 1,8                        | 959994-902                    |                              |                                         |                          |                              |                             |                      |
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 100               | 1,8                        | 959964-902                    | 959964-906                   | 959964-912                              | 959964-918               | 928975-902                   |                             | 728975-902           |
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 75                | 1,8                        | 959951-902                    |                              |                                         |                          |                              |                             |                      |
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 50                | 1,8                        | 959941-902                    | 959941-906                   | 959941-912                              | 959941-918               | 927975-902                   | 927975-906                  | 727975-902           |
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 30                | 1,8                        | 959931-902                    | 959931-906                   | 959931-912                              | 959931-918               | 924975-902                   | 924975-906                  | 724975-902           |
|                                                                                     | Rapid Resolution<br>HT, 600 bar          | 4,6 x 20                | 1,8                        |                               |                              |                                         |                          | 926975-902                   | 926975-906                  | 726975-902           |
|    | UHPLC-Vorsäule,<br>600 bar,<br>3 St./VE  | 4,6 x 5                 | 1,8                        | 820750-901                    |                              |                                         |                          | 820750-903                   |                             |                      |
|                                                                                     | Solvent Saver<br>HT, 600 bar             | 3,0 x 100               | 1,8                        | 959964-302                    | 959964-306                   | 959964-312                              |                          | 928975-302                   |                             | 728975-302           |
|                                                                                     | Solvent Saver<br>HT, 600 bar             | 3,0 x 50                | 1,8                        | 959941-302                    | 959941-306                   | 959941-312                              |                          | 927975-302                   | 927975-306                  | 727975-302           |
|                                                                                     | Solvent Saver<br>HT, 600 bar             | 3,0 x 30                | 1,8                        |                               |                              |                                         |                          | 924975-302                   | 924975-306                  | 724975-302           |
|                                                                                     | Solvent Saver<br>HT, 600 bar             | 3,0 x 20                | 1,8                        |                               |                              |                                         |                          | 926975-302                   | 926975-306                  | 726975-302           |
|  | UHPLC-Vorsäule,<br>1200 bar,<br>3 St./VE | 3,0 x 5                 | 1,8                        | 823750-901                    |                              |                                         |                          | 823750-903                   |                             |                      |
|                                                                                     | Narrowbore<br>RRHT, 600 bar              | 2,1 x 150               | 1,8                        | 959794-902                    |                              |                                         |                          |                              |                             |                      |
|                                                                                     | Narrowbore<br>RRHT, 600 bar              | 2,1 x 100               | 1,8                        | 959764-902                    | 959764-906                   | 959764-912                              | 959764-918               | 928700-902                   | 928700-906                  | 728700-902           |
|                                                                                     | Narrowbore<br>RRHT, 600 bar              | 2,1 x 50                | 1,8                        | 959741-902                    | 959741-906                   | 959741-912                              | 959741-918               | 927700-902                   | 927700-906                  | 727700-902           |
|                                                                                     | Narrowbore<br>RRHT, 600 bar              | 2,1 x 30                | 1,8                        | 959731-902                    | 959731-906                   | 959731-912                              |                          | 924700-902                   | 924700-906                  | 724700-902           |
|                                                                                     | Narrowbore<br>RRHT, 600 bar              | 2,1 x 20                | 1,8                        |                               |                              |                                         |                          | 926700-902                   | 926700-906                  | 726700-902           |
|  | UHPLC-Vorsäule,<br>1200 bar,<br>3 St./VE | 2,1 x 5                 | 1,8                        | 821725-901                    |                              |                                         |                          | 821725-903                   |                             |                      |