

LC/MS-Kalibriermischung und Quellen-Kompatibilitätsmatrix

Quelle	Serie 6100 Single Quadrupol LC/MS*	6140A/6150B Single Quadrupol LC/MS	Serie 6200 TOF LC/MS	Serie 6300 Ion Trap LC/MS	Serie 6400 Triple Quadrupol LC/MS	Serie 6500 Accurate-Mass Q-TOF LC/MS
ESI	G2421-60001	G1969-85000	G1969-85000	G2431A	G1969-85000	G1969-85000
APCI	G2432A	G1969-85010	G1969-85010	G2432A	G1969-85010 ²	G1969-85010 ²
APPI	G2432A	G2432A	G1969-85010	G2432A	G2432A ²	G1969-85010 ²
MMI	G2432A	G1969-85000	G1969-85020	G2432A	G1969-85020	G1969-85020
NanoESI			G1969-85000 ⁴	G2431A		G1969-85000 ⁴
HPLC-Chip Cube		G1969-85000 ²	G1969-85000 ⁴	G2431A ¹	G1969-85000 ³	G1969-85000 ⁴

*G6110A, G6120A/B, G6130A/B

¹1:5-Verdünnung empfohlen²Kein Autotune³Nur positives ESI-Tuning⁴Nur Kalibrierung