

20 % Sebaconitril/2 % H₃PO₄

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,61 m (2 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80015	G3591-81015	G3591-82015
9,14 m (30 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80014	G3591-81014	G3591-82014

Silicagel

Länge	AD (Zoll)	ID (mm)	Mesh	Edelstahl	UltiMetal	Nickel	
0,61 m (2 ft)		1/8	2	60/80	G3591-80141	G3591-81141	G3591-82141
1,22 m (4 ft)		1/8	2	60/80	G3591-80142	G3591-81142	G3591-82142
1,83 m (6 ft)		1/8	2	60/80	G3591-80108		

0,1 % SP-1000

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
2,13 m (7 ft)	1/8	2	Carbopak C	80/100	G3591-80063	G3591-81063	G3591-82063

25 % SP-2100

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,5 m (1,64 ft)	1/16	1	Chromosorb PAW	80/100	G3591-80007	G3591-81007	G3591-82007
1,75 m (5,7 ft)	1/16	1	Chromosorb PAW	80/100	G3591-80008	G3591-81008	G3591-82008
4,57 m (15 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80068	G3591-81068	G3591-82068

20 % TCEP

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,56 m (1,84 ft)	1/16	1	Chromosorb PAW	80/100	G3591-80006	G3591-81006	G3591-82006
1,52 m (5 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80094	G3591-81094	G3591-82094
4,57 m (15 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80049	G3591-81049	G3591-82049

10 % UC W982

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,46 m (1,5 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80034	G3591-81034	G3591-82034
0,61 m (2 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80040	G3591-81040	G3591-82040

12 % UC W982

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,61 m (2 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80000	G3591-81000	G3591-82000