

10% Carbowax 20M (G16, G\$1)

| Length | OD (in) | ID (mm) | Support | Mesh | Stainless Steel |
|---------------|---------|---------|----------------|--------|-----------------|
| 6.56 ft (2 m) | 1/8 | 2 | Chromosorb WHP | 80/100 | G3591-70016 |

10% Carbowax 20M (G16, G\$1) + 2% KOH

| Length | OD (in) | ID (mm) | Support | Mesh | Stainless Steel |
|-----------------|---------|---------|----------------|--------|-----------------|
| 5.91 ft (1.8 m) | 1/8 | 2 | Chromosorb WHP | 80/100 | G3591-70012 |

20% Carbowax 20M (G16, G\$1)

| Length | OD (in) | ID (mm) | Support | Mesh | UltiMetal | Stainless Steel | Nickel |
|---------------|---------|---------|----------------|---------|-------------|-----------------|-------------|
| 9.84 ft (3 m) | 1/8 | 2 | Chromosorb WHP | 100/120 | G3591-81099 | G3591-80099 | G3591-82099 |

7% Carbowax M + 3% Polyphenoether 6 ring + 2% KOH

| Length | OD (in) | ID (mm) | Support | Mesh | UltiMetal | Nickel |
|---------------|---------|---------|----------------|--------|-------------|-------------|
| 4 ft (1.22 m) | 1/8 | 2 | Chromosorb WAW | 80/100 | G3591-81050 | G3591-82050 |

Carboxen-1000

| Length | OD (in) | ID (mm) | Mesh | UltiMetal | Stainless Steel |
|----------------|---------|---------|-------|-------------|-----------------|
| 10 ft (3.05 m) | 1/8 | 2 | 60/80 | G3591-81055 | G3591-80055 |

Chromosorb 101

| Length | OD (in) | ID (mm) | Mesh | UltiMetal | Stainless Steel |
|---------------|---------|---------|--------|-------------|-----------------|
| 6 ft (1.83 m) | 1/8 | 2 | 80/100 | G3591-81021 | G3591-80021 |

Chromosorb 102

| Length | OD (in) | ID (mm) | Mesh | UltiMetal | Stainless Steel | Nickel |
|---------------|---------|---------|--------|-------------|-----------------|-------------|
| 2 ft (0.61 m) | 1/8 | 2 | 80/100 | G3591-81139 | G3591-80139 | G3591-82139 |