

**Porapak Q**

Länge	AD (Zoll)	ID (mm)	Mesh	Edelstahl	UltiMetal	Nickel
0,91 m (3 ft)	1/8	2	80/100	G3591-80135	G3591-81135	G3591-82135
1,83 m (6 ft)	1/8	2	60/80	G3591-80136	G3591-81136	G3591-82136
1,83 m (6 ft)	1/8	2	80/100	G3591-80005	G3591-81005	G3591-82005
1,83 m (6 ft)	1/8	2	80/100	G3591-80013	G3591-81013	G3591-82013
2,44 m (8 ft)	1/8	2	60/80	G3591-80137	G3591-81137	G3591-82137
2,5 m (8,2 ft)	1/8	2	80/100	G3591-80083	G3591-81083	G3591-82083
2,74 m (9 ft)	1/8	2	80/100	G3591-80016	G3591-81016	G3591-82016
3,96 m (13 ft)	1/8	2	80/100	G3591-80053	G3591-81053	G3591-82053
4,57 m (15 ft)	1/8	2	80/100	G3591-80066	G3591-81066	G3591-82066
7,62 m (25 ft)	1/8	2	100/120	G3591-80052	G3591-81052	G3591-82052
9,14 m (30 ft)	1/16	1	80/100	G3591-80096	G3591-81096	G3591-82096
1 m (3,28 ft.)	1/8	2	80/100	G3591-70014		
1,8 m (5,91 ft.)	1/8	2	80/100	G3591-70010		
3 m (9,84 ft.)	1/8	2	80/100	G3591-70009		
2 m (6,56 ft.)	1/8	2	80/100	G3591-70001		
1,5 m (4,92 ft.)	1/8	2	50/80	G3591-70018		

**Porapak QS**

Länge	AD (Zoll)	ID (mm)	Mesh	Edelstahl	UltiMetal	Nickel
2 m (6,56 ft)	1/8	2	80/100	G3591-80157		
2,44 m (8 ft)	1/8	2	80/100	G3591-80051	G3591-81051	G3591-82051
1,8 m (5,91 ft.)	1/8	2	80/100	G3591-70010		
3 m (9,84 ft.)	1/8	2	80/100	G3591-70009		

**Porapak R**

Länge	AD (Zoll)	ID (mm)	Mesh	Edelstahl	UltiMetal	Nickel
1,83 m (6 ft)	1/8	2	60/80	G3591-80106	G3591-81106	G3591-82106

**Porapak T**

Länge	AD (Zoll)	ID (mm)	Mesh	Edelstahl	UltiMetal	Nickel
0,46 m (1,5 ft)	1/8	2	80/100	G3591-80138	G3591-81138	G3591-82138
2 m (6,56 ft)	1/8	2	80/100	G3591-80120	G3591-81120	G3591-82120

**20 % Sebaconitril**

Länge	AD (Zoll)	ID (mm)	Träger	Mesh	Edelstahl	UltiMetal	Nickel
0,61 m (2 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80029	G3591-81029	G3591-82029
6 m (19,7 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80071	G3591-81071	G3591-82071
9,14 m (30 ft)	1/8	2	Chromosorb PAW	80/100	G3591-80026	G3591-81026	G3591-82026