

Präparierbesteck, 8-teilig

1 Schere spitz	A=110
1 Pinzette anatomisch	A=145
1 Skalpell 35 mm Klinge	A=150
1 Skalpell 45 mm Klinge	A=160
1 Verbandschere	A=140
1 Pinzette spitz	A=120
1 Präpariernadel	A=140
1 Lanzette	A=140

Dissecting set, 8 pieces

1 Scissors, sharp	A=110
1 Forceps, blunt	A=145
1 Scalpel 35 mm blade	A=150
1 Scalpel 45 mm blade	A=160
1 Dressing scissors	A=140
1 Forceps, sharp	A=120
1 Dissection needle	A=140
1 Lancet	A=140



12017

Material	Artikel Nr.
Rostfrei	12017



33999

Labor-Set, 4-teilig

Ideal für den Einstieg in die Laborarbeit.

4 Instrumente im praktischen Labor-Set bestehend aus:

- 1 Verbandschere (4100)
- 1 Mikrodoppelspatel (3017)
- 1 Pinzette, spitz (1012)
- 1 Mikroschaufel, halbrund (3303)

Laboratory-Set, 4 pieces

Ideal for the first works in the laboratory.

4 instruments in a practical Laboratory-Set consisting in:

- 1 Dressing scissor (4100)
- 1 Micro spatula, double (3017)
- 1 Forceps, sharp (1012)
- 1 Micro scoop, half rounded (3303)

Material	A	B	Artikel Nr.
18/10-Stahl	170	120	33999

Präpariernadeln mit Edelstahlgriff

Typ 1: gerade, Typ 2: gebogen, Typ 3: lanzettenförmig

Dissecting needle with stainless steel handle

Type 1: straight, Type 2: bent, Type 3: lancet

Material	A	Typ	Artikel Nr.
18/10-Stahl	140	1	12030
18/10-Stahl	140	2	12035
18/10-Stahl	140	3	12020



Typ 1

12030

Typ 2

12035

Typ 3

12020