


Doppelspatel, Typ Chattaway Double spatulas, type Chattaway

| Material | A | a | B | D | Artikel Nr. |
|-------------|-----|----|----|-----|-------------|
| 18/10-STAHL | 125 | 40 | 7 | 3 | 3141 |
| 18/10-STAHL | 150 | 40 | 7 | 3,5 | 3142 |
| 18/10-STAHL | 180 | 50 | 8 | 3,5 | 3143 |
| 18/10-STAHL | 200 | 60 | 9 | 4 | 3144 |
| 18/10-STAHL | 235 | 70 | 10 | 4 | 3145 |



Doppelspatel, 1 Seite konisch Double spatulas, 1 side tapered

| Material | A | a | B | D | Artikel Nr. |
|--------------|-----|----|---|-----|-------------|
| 18/10-Stahl | 210 | 50 | 7 | 2,5 | 3170 |
| 18/10-Stahl | 160 | 50 | 9 | 3,5 | 3175 |
| PTFE Coating | 210 | 50 | 7 | 2,5 | 3730 |

PTFE = Polytetrafluorethylen/polytetrafluorethylene

Typ 1

3070

Typ 2

3071

Typ 3

3072

Zementspatel Cement spatulas

| Material | A | a | D | Typ | Artikel Nr. |
|----------|-----|----|---|-----|-------------|
| Rostfrei | 150 | 35 | 5 | 1 | 3070 |
| Rostfrei | 150 | 35 | 5 | 2 | 3071 |
| Rostfrei | 150 | 35 | 5 | 3 | 3072 |



Mörser Doppelspatel

Massive Ausführung

Mortar spatulas

Heavy type

| Material | A | a | B | D | Artikel Nr. |
|-------------|-----|----|----|----|-------------|
| 18/10-Stahl | 150 | 40 | 20 | 8 | 3260 |
| 18/10-Stahl | 200 | 40 | 20 | 8 | 3261 |
| 18/10-Stahl | 250 | 60 | 25 | 10 | 3262 |
| 18/10-Stahl | 300 | 60 | 25 | 10 | 3263 |
| 18/10-Stahl | 400 | 60 | 25 | 10 | 3264 |