




## Labor-Scheren Laboratory scissors



Typ 1: rund-rund  
Typ 2: spitz-rund  
Typ 3: spitz-spitz

Typ 1: round-round  
Typ 2: sharp-round  
Typ 3: sharp-sharp



### Laborscheren mit Kunststoffgriff Laboratory scissors with plastic handle

Material	A	a	Typ	Artikel Nr.
Rostfrei	130	50	2	4020
Rostfrei	150	70	2	4021
Rostfrei	200	100	2	4022
Rostfrei	250	120	2	4023



4020



4070



4071



4080

### Mikroskopierscheren

Schere gerade (4070, 4080) und gebogen (4071)

### Scissors microscopy

Scissor straight (4070, 4080) and bent (4071)

Material	A	a	Typ	Artikel Nr.
Rostfrei	100	25	3	4070
Rostfrei	100	30	3	4071
Rostfrei	115	35	3	4080