


5300

Doppelmuffe, 18/10-Stahl, DIN 12895
Bosshhead, 18/10 stainless steel, DIN 12895

| Material | d | M | Typ | Winkel ↖ | Artikel Nr. |
|-------------|------|---|-----|----------|-------------|
| 18/10-Stahl | 16,5 | 6 | 2 | 90 | 5300 |
| 18/10-Stahl | 16,5 | 6 | 3 | 90 | 5310 |
| 18/10-Stahl | 16,5 | 8 | 4 | 90 | 5305 |



5360

Doppelmuffe, T-Guss, DIN 12895

Verchromt

Bosshhead, malleable cast iron, DIN 12895

Chromium plated

| Material | d | M | Typ | Winkel ↖ | Artikel Nr. |
|------------------|------|----|-----|----------|-------------|
| T-Guss verchromt | 16,5 | 8 | 1 | 90 | 5360 |
| T-Guss verchromt | 16,5 | 8 | 4 | 90 | 5365 |
| T-Guss verchromt | 20,0 | 8 | 1 | 90 | 5359 |
| T-Guss verchromt | 36,0 | 10 | 1 | 90 | 5369 |

5332



Doppelmuffe, Aluminium, DIN 12895
Bosshhead, aluminium, DIN 12895

| Material | d | M | Typ | Winkel ↖ | Artikel Nr. |
|----------|------|---|-----|----------|-------------|
| Alu | 16,5 | 8 | 1 | 90 | 5332 |
| Alu | 16,5 | 8 | 3 | 90 | 5342 |
| Alu | 16,5 | 8 | 4 | 90 | 5337 |

