




5460

Mikroklemme, T-Guss

T-Guss grün lackiert, Finger mit Kork.

Clamp micro, malleable cast iron

Malleable cast iron green varnished, finger with cork.



Material	A	D	d	Artikel Nr.
T-Guss	150	12	15	5460




5475

Mikroklemme, Stahl verzinkt

Finger mit Silikon

Clamp micro, steel zincked

Finger with silicone



Material	A	D	d	Artikel Nr.
Stahl verzinkt	150	12	0-20	5475

Universalklemme DIN 12894, Aluminium

Finger mit Kork

Universal clamp DIN 12894, aluminium

Finger with cork



5444

Material	A	D	d	Artikel Nr.
Alu	180	12	0-80	5444




5425

Universalklemme DIN 12894, 18/10-Stahl

Finger mit Kork

Universal clamp DIN 12894, 18/10 stainless steel

Finger with cork



Material	A	D	d	Artikel Nr.
18/10-Stahl	150	12	0-80	5425