

## BUNSEN BURNER



Lab Bunsen burner for butane gas with suitable cartridges and for connecting to a gas pipe.

ITEM NO.	DESCRIPTION	MANUFACTURER	QUAN-TITY	UNIT
<b>HE1395</b>	Labogaz Bunsen burner for butane gas, without cartridge, height incl. cartridge 24 cm	Hecht	1	piece
<b>HE1395/1</b>	Cartridge with 450 g (806 ml) Butane gas CV 470 for Labogaz burner 1395. With safety valve (cartridge can be separated from the burner at any time without risk, when the unit is not operating)	Hecht	1	piece
<b>HE1396</b>	LABOGAZ Bunsen burner, without cartridge. Mounting system for lancets. A plastic connection prevents opening of the metal brackets and therefore a cartridge change without removal of the device	Hecht	1	piece
<b>HE1396/1</b>	Butane gas cartridge for Labogaz burner 1396, 190 g	Hecht	1	piece
<b>BO7360</b>	All gas Bunsen burner with air regulation, connection to gas line, with tilted tap, max 1300°C, air regulation, ignition flame, KW 0,5-3	Bochem	1	piece