

2017 Product List

| REFERENCE | PHASE | LENGTH (mm) | ID (mm) | PARTICLE SIZE (µ) | TYPE |
|-----------|--------------------|----------------|------------|----------------------|------------|
| 52515 | CHIRALPAK® ZWIX(-) | 250 | 4 | 3 | Analytical |
| 52584 | CHIRALPAK® ZWIX(-) | 150 | 3 | 3 | Analytical |
| 52585 | CHIRALPAK® ZWIX(-) | 250 | 3 | 3 | Analytical |

CHIRALPAK® QD-AX & QN-AX

| | | | | | |
|-------|----------------------|-----|-----|---|------------------------|
| 31311 | CHIRALPAK® QD-AX | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 31324 | CHIRALPAK® QD-AX | 150 | 4.6 | 5 | Analytical |
| 31394 | CHIRALPAK® QD-AX | 150 | 2.1 | 5 | Analytical |
| 31344 | CHIRALPAK® QD-AX | 150 | 20 | 5 | Semi-prep |
| 31444 | CHIRALPAK® QD-AX SFC | 150 | 20 | 5 | Semi-prep |
| 32311 | CHIRALPAK® QN-AX | 10 | 4 | 5 | Guard Cartridges (x 3) |
| 32324 | CHIRALPAK® QN-AX | 150 | 4.6 | 5 | Analytical |
| 32394 | CHIRALPAK® QN-AX | 150 | 2.1 | 5 | Analytical |
| 32344 | CHIRALPAK® QN-AX | 150 | 20 | 5 | Semi-prep |
| 32444 | CHIRALPAK® QN-AX SFC | 150 | 20 | 5 | Semi-prep |

CHIRALPAK® MA(+)

| | | | | | |
|-------|------------------|----|-----|---|------------|
| 21822 | CHIRALPAK® MA(+) | 50 | 4.6 | 3 | Analytical |
|-------|------------------|----|-----|---|------------|

CHIRALPAK® WH

| | | | | | |
|-------|---------------|-----|-----|----|------------|
| 25622 | CHIRALPAK® WH | 50 | 4.6 | 10 | Guard |
| 25625 | CHIRALPAK® WH | 250 | 4.6 | 10 | Analytical |

SFC ACHIRAL COLUMNS

| | | | | | |
|-------|-------------|-----|-----|---|------------|
| AA424 | DCpak SFC-A | 150 | 4.6 | 5 | Analytical |
| AA425 | DCpak SFC-A | 250 | 4.6 | 5 | Analytical |
| AA435 | DCpak SFC-A | 250 | 10 | 5 | Analytical |
| AA445 | DCpak SFC-A | 250 | 20 | 5 | Analytical |
| AA494 | DCpak SFC-A | 150 | 2.1 | 5 | Analytical |
| AA495 | DCpak SFC-A | 250 | 2.1 | 5 | Analytical |
| AB424 | DCpak SFC-B | 150 | 4.6 | 5 | Analytical |
| AB425 | DCpak SFC-B | 250 | 4.6 | 5 | Analytical |
| AB435 | DCpak SFC-B | 250 | 10 | 5 | Analytical |
| AB445 | DCpak SFC-B | 250 | 20 | 5 | Analytical |
| AB494 | DCpak SFC-B | 150 | 2.1 | 5 | Analytical |
| AB495 | DCpak SFC-B | 250 | 2.1 | 5 | Analytical |

COLUMN ACCESSORIES

| | | | | | |
|------------|--|-----|-----|---|-----------------------------|
| 00011 | GUARD CARTRIDGE HOLDER† | - | - | - | Holder for guard cartridges |
| 00024 | Column Jacket | 150 | 4.6 | - | Hardware |
| 00025 | Column Jacket | 250 | 4.6 | - | Hardware |
| 000D1-DAIC | Guard Column Coupler | - | - | - | Hardware |
| FK002 | Guard Filter Kit (1 holder, 3 glass filters, 2 fittings and 1 tubing) | - | - | - | Guard Filter |
| FR002 | Replacement Guard Filters (x 6) | - | - | - | Guard Filters |

†For all guard cartridges