

LABOMECA® 3- and 4-Finger Clamps, BOCHEM

Three-finger clamps

	Material	Shaft	Clamp width mm	Rod Ø mm	Rod L mm	Order No.
①	Aluminium alloy	flexible	0 – 30	8	220	31 05 26007
②	Aluminium alloy/PVC ¹⁾		12 – 100	12	140	31 05 26003

Four-finger clamps

	Material	Clamp width mm	Rod Ø mm	Rod L mm	Order No.
③	18/10 steel/silicone ¹⁾	35 – 100	12	140	31 05 26004

¹⁾ Span width independently adjustable by 2 regulating screws.

