

## LABOMECA® Clamps, BOCHEM

### Clamps, round jaws, plastic coating

	Material	Clamp width mm	Rod Ø mm	Rod L mm	Order No.
①	LABORAL alloy/PVC	8 – 30	8	125	31 05 27005
②	18/10 steel/silicone	8 – 30	8	125	31 05 27006
①	LABORAL alloy/PVC	15 – 65	12	140	31 05 27002
②	18/10 steel/silicone	20 – 75	12	140	31 05 27003
②	LABORAL alloy/PVC	22 – 70	12	140	31 05 27004

### Clamp, flexible shaft, plastic coating

	Material	Clamp width mm	Rod Ø mm	Rod L mm	Order No.
③	Aluminium alloy	15 – 65	12	225	31 05 27001

