

## Reusable Bottletop Filters, NALGENE®, IDL

### ① Polysulfone (PSF) bottletop filters with silicone O-Ring

Can be autoclaved repeatedly, intended for use with threaded 33 or 45 mm dia. media bottles. Suitable for pressure and vacuum filtration up to 0.6 bar. Non-breakable, non-toxic. Includes seal, cap and tubing adapter made of polypropylene. 2 filter upper sizes available (250 or 500 ml), which accept 47-mm or 50-mm membranes.

Volume top part ml	Thread mm	Order No.
500	33	26 87 32133
500	45	26 87 32145

### ② Polysulfone (PSF) filter holder with filtrate bottle

You will no longer need to use stainless-steel vessels or breakable glass Büchner flasks. Graded filtrate bottle. Two connectors for connecting the tubes to negative pressure or syringe filter. A cotton plug can be inserted into the hose adapter for sterile aeration. Tubes of an inner diameter of 6 to 8 mm can be connected to the adapter. Closed with a cap, the sterile filtrate can be stored in the bottle for greater user comfort. Autoclavable, graded, patented.

Volume top part ml	Volume filtrate bottle ml	Overall width incl. side-arms mm	Filtrate bottle Ø mm	Height incl. cover ports mm	Order No.
250	250	110	83	180	10 87 30400
500	500	134	117	230	10 87 30405
500	1000	134	114	293	10 87 30410

### ③ Polysulfone (PSF) filter holder with funnel

Suitable for use on filter bottles or suction devices together with an appropriate rubber stopper with bore hole. Supplied with an easy-to-use, non-twisting NALGENE vacuum seal that may be used instead of a rubber stopper no. 7 or no. 8 with 15 mm bore hole. Autoclavable, patented.

Volume top part ml	Ø mm	Height incl. cover ports mm	Order No.
250	94	190	10 87 31400
500	120	214	10 87 31405

