

Bulk sorbents for preparation of TLC plates

Ordering information – Aluminium oxides for TLC and PLC (particle size 5 – 40 µm)

| Product | Ordering No. | pH of 10 % aqueous suspension | Package | Contents of one package | Method |
|---|--------------|-------------------------------|---------|-------------------------|--------|
| Aluminium oxide 60 G neutral | 1.01090.2500 | 7.5 | Plastic | 2.5 kg | TLC |
| | 1.01090.9025 | 7.5 | Plastic | 25 kg | TLC |
| Aluminium oxide 60 G F ₂₅₄ neutral | 1.01092.0500 | 7.5 | Plastic | 500 g | TLC |

F₂₅₄: Green fluorescent indicator

Ordering information – Other sorbents for TLC

| Product | Ordering No. | Particle size | Package | Contents of one package |
|----------------------------|--------------|---------------|---------|-------------------------|
| Cellulose microcrystalline | 1.02330.0500 | <20 µm | Plastic | 500 g |

