

Ordering information – Purospher® STAR RP-18 endcapped (3 µm), stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Purospher® STAR RP-18 endcapped, cartridge set (1 LiChroCART® 30-2 and 1 manu-CART® 30 mm)	1.50237.0001	3 µm	30 mm	2 mm	1 set
Purospher® STAR RP-18 endcapped	1.50238.0001	3 µm	30 mm	2 mm	3 pieces
Purospher® STAR RP-18 endcapped, cartridge set (1 LiChroCART® 55-2 and 1 manu-CART® 55 mm)	1.50240.0001	3 µm	55 mm	2 mm	1 set
Purospher® STAR RP-18 endcapped, cartridge set	1.50241.0001	3 µm	55 mm	2 mm	3 pieces
Purospher® STAR RP-18 endcapped, custom packed	1.50190.7184	3 µm	125 mm	2 mm	1 piece
Purospher® STAR RP-18 endcapped, custom packed	1.50175.7184	3 µm	125 mm	3 mm	1 piece
Purospher® STAR RP-18 endcapped, custom packed	1.50177.7184	3 µm	250 mm	3 mm	1 piece
Purospher® STAR RP-18 endcapped, cartridge set (1 LiChroCART® 30-4 and 1 manu-CART® 30 mm)	1.50239.0001	3 µm	30 mm	4 mm	1 set
Purospher® STAR RP-18 endcapped	1.50225.0001	3 µm	30 mm	4 mm	3 pieces
Purospher® STAR RP-18 endcapped, cartridge set (1 LiChroCART® 55-4 and 1 manu-CART® 55 mm)	1.50242.0001	3 µm	55 mm	4 mm	1 set
Purospher® STAR RP-18 endcapped	1.50231.0001	3 µm	55 mm	4 mm	3 pieces
Purospher® STAR RP-18 endcapped	1.51460.0001	3 µm	75 mm	4 mm	1 piece
Purospher® STAR RP-18 endcapped, custom packed	1.50174.7184	3 µm	250 mm	4 mm	1 piece
Purospher® STAR RP-18 endcapped, custom packed	1.51448.7184	3 µm	100 mm	4.6 mm	1 piece

The LiChroCART® columns (75, 125, 150 and 250 mm length) in the list above (2, 3, 4 and 4.6 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. LiChroCART® columns 250–10 mm require part number 1.51419.0001 manu-CART® 10. The short LiChroCART® columns (30 and 55 mm length) can be ordered as a set including the corresponding cartridge holder and one cartridge, or as a pack of 3 cartridges without cartridge holder. Additional dimensions available as customized packings see page 366. The separate part numbers for the cartridge are as follows: 1.50227.0001 LiChroCART® cartridge holder for 30 mm cartridge and 1.50226.0001 LiChroCART® cartridge holder for 55 mm cartridge.