

## Ordering information – LiChrospher® 100 RP-18, stainless steel columns Hibar®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
LiChrospher® 100 RP-18	1.00423.7079	5 µm	250 mm	3 mm	1 piece
LiChrospher® 100 RP-18	1.50477.0001	5 µm	125 mm	4 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50477.1003	5 µm	125 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50377.0001	5 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50377.1003	5 µm	250 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50545.0001	5 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50545.1003	5 µm	100 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.50546.0001	5 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50546.1003	5 µm	150 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.50547.0001	5 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC columns from 3 batches)	1.50547.1003	5 µm	250 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.50183.7079	5 µm	250 mm	10 mm	1 piece
LiChrospher® 100 RP-18	1.50181.7080	7 µm	125 mm	4 mm	1 piece
LiChrospher® 100 RP-18	1.50012.7080	7 µm	125 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.50009.7080	7 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.00424.7080	7 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.50013.7081	10 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.00424.7081	10 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.50183.7081	10 µm	250 mm	10 mm	1 piece

The Hibar® columns are complete with endfittings. When using a guard column with a Hibar® column, we recommend part number 1.51487.0001 guard column cartridge holder for 4–4 mm guard column cartridges LiChroCART®. Additional dimensions available as customized packings see page 366.

## Ordering information – LiChrospher® 100 RP-18 endcapped, stainless steel columns Hibar®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
LiChrospher® 100 RP-18 endcapped	1.51907.0001	5 µm	125 mm	2 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50182.7085	5 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50548.0001	5 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.51906.0001	5 µm	125 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50549.0001	5 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50550.0001	5 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50183.7085	5 µm	250 mm	10 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50182.7084	10 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50012.7084	10 µm	125 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50009.7084	10 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.00424.7084	10 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 endcapped	1.50183.7084	10 µm	250 mm	10 mm	1 piece