

# LiChrospher® 100 RP-18 and RP-18 endcapped

## Ordering information – LiChrospher® 100 RP-18, stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
LiChrospher® 100 RP-18	1.50195.7079	5 µm	125 mm	2 mm	1 piece
LiChrospher® 100 RP-18	1.50190.7079	5 µm	250 mm	2 mm	1 piece
LiChrospher® 100 RP-18	1.50233.7079	5 µm	30 mm	3 mm	1 piece
LiChrospher® 100 RP-18	1.50159.0001	5 µm	125 mm	3 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50159.1003	5 µm	125 mm	3 mm	3 pieces
LiChrospher® 100 RP-18	1.50154.0001	5 µm	250 mm	3 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50154.1003	5 µm	250 mm	3 mm	3 pieces
LiChrospher® 100 RP-18	1.50957.0001	5 µm	4 mm	4 mm	10 pieces
LiChrospher® 100 RP-18	1.50931.0001	5 µm	25 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50302.7079	5 µm	30 mm	4 mm	1 piece
LiChrospher® 100 RP-18	1.50228.7079	5 µm	55 mm	4 mm	1 piece
LiChrospher® 100 RP-18	1.50987.0001	5 µm	75 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50823.0001	5 µm	125 mm	4 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50823.1003	5 µm	125 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50943.0001	5 µm	125 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50833.0001	5 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50833.1003	5 µm	250 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50983.0001	5 µm	250 mm	4 mm	3 pieces
LiChrospher® 100 RP-18	1.50600.0001	5 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50600.1003	5 µm	100 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.50601.0001	5 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50601.1003	5 µm	150 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.50602.0001	5 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18 (3 HPLC cartridges from 3 batches)	1.50602.1003	5 µm	250 mm	4.6 mm	3 pieces
LiChrospher® 100 RP-18	1.51442.7080	7 µm	125 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.51432.7080	7 µm	150 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.51431.7080	7 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.50843.0001	10 µm	250 mm	4 mm	1 piece
LiChrospher® 100 RP-18	1.51448.7081	10 µm	100 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.51431.7081	10 µm	250 mm	4.6 mm	1 piece
LiChrospher® 100 RP-18	1.50853.0001	10 µm	250 mm	10 mm	1 piece

The LiChroCART® columns (75, 100, 125, 150 and 250 mm length) in the list above (2, 3, 4 and 4.6 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. LiChroCART® columns 250–10 mm require part number 1.51419.0001 manu-CART® 10. Additional dimensions and validation kit available as customized packings see page 366.