

Purospher® STAR Si and Purospher® STAR NH₂

Ordering information – Purospher® STAR Si, stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Purospher® STAR Si, custom packed	1.50175.7174	3 µm	125 mm	3 mm	1 piece
Purospher® STAR Si, custom packed	1.50177.7174	3 µm	150 mm	3 mm	1 piece
Purospher® STAR Si, custom packed	1.50190.7175	5 µm	150 mm	2 mm	1 piece
Purospher® STAR Si, custom packed	1.50175.7175	5 µm	125 mm	3 mm	1 piece
Purospher® STAR Si	1.50249.0001	5 µm	4 mm	4 mm	10 pieces
Purospher® STAR Si	1.50268.0001	5 µm	125 mm	4 mm	1 piece
Purospher® STAR Si	1.50269.0001	5 µm	250 mm	4 mm	1 piece
Purospher® STAR Si	1.50356.0001	5 µm	150 mm	4.6 mm	1 piece
Purospher® STAR Si	1.50357.0001	5 µm	250 mm	4.6 mm	1 piece

The LiChroCART® columns (75, 125, 150 and 250 mm length) in the list above (4 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. Additional dimensions available as customized packings see page 366.

Ordering information – Purospher® STAR Si, stainless steel columns Hibar® RT

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Purospher® STAR Si, custom packed	1.50181.7174	3 µm	125 mm	4 mm	1 piece
Purospher® STAR Si, custom packed	1.50182.7174	3 µm	250 mm	4 mm	1 piece
Purospher® STAR Si, custom packed	1.50012.7174	3 µm	125 mm	4.6 mm	1 piece
Purospher® STAR Si, custom packed	1.50009.7174	3 µm	150 mm	4.6 mm	1 piece
Purospher® STAR Si, custom packed	1.50009.7175	5 µm	150 mm	4.6 mm	1 piece

The Hibar® RT columns are complete with endfittings. When using a guard column with a Hibar® RT column, we recommend part number 1.51487.0001 guard column cartridge holder for 4–4 mm guard column cartridges LiChroCART®. Additional dimensions available as customized packings see page 366.

Ordering information – Purospher® STAR NH₂, stainless steel cartridges LiChroCART®

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Purospher® STAR NH ₂	1.50190.7177	5 µm	250 mm	2 mm	1 piece
Purospher® STAR NH ₂	1.50177.7177	5 µm	250 mm	3 mm	1 piece
Purospher® STAR NH ₂	1.50267.0001	5 µm	4 mm	4 mm	10 pieces
Purospher® STAR NH ₂	1.50244.0001	5 µm	125 mm	4 mm	1 piece
Purospher® STAR NH ₂	1.50245.0001	5 µm	250 mm	4 mm	1 piece
Purospher® STAR NH ₂	1.50247.0001	5 µm	150 mm	4.6 mm	1 piece
Purospher® STAR NH ₂	1.50248.0001	5 µm	250 mm	4.6 mm	1 piece

The LiChroCART® columns (75, 125, 150 and 250 mm length) in the list above (4 mm i.d.) require part number 1.51486.0001 manu-CART® cartridge column holder, which can be used to hold one cartridge column with or without a 4–4 mm guard column. Additional dimensions available as customized packings see page 366.

Ordering information – Purospher® STAR NH₂, stainless steel columns Hibar® RT

Product	Ordering No.	Particle size	Dimension length	Dimension i.d.	Contents of one package
Purospher® STAR NH ₂ , custom packed	1.50182.7177	5 µm	250 mm	4 mm	1 piece
Purospher® STAR NH ₂ , custom packed	1.50009.7177	5 µm	150 mm	4.6 mm	1 piece

The Hibar® RT columns are complete with endfittings. When using a guard column with a Hibar® RT column, we recommend part number 1.51487.0001 guard column cartridge holder for 4–4 mm guard column cartridges LiChroCART®. Additional dimensions available as customized packings see page 366.