



Performance

Jumbosep™ Device Reduces Processing Time Over Other Devices



1 mg/mL BSA solution was processed in each of the above devices until a 15-fold concentration was achieved.

Concentration Time



Concentrate dilute protein samples in less than 30 minutes with 10, 30, and 100K Jumbosep devices.

Ordering Information

The generic starter kit includes four holders, cups, and caps. Membrane inserts sold separately. The starter kits include four holders, cups, caps, and membrane inserts.

Jumbosep Centrifugal Device Starter Kits

Part Number	Description	Pkg
FD000K65	Generic starter kit, (no membrane inserts)	4/pkg
FD003K65	3K starter kit, gray	4/pkg
FD010K65	10K starter kit, blue	4/pkg
FD030K65	30K starter kit, red	4/pkg
FD100K65	100K starter kit, clear	4/pkg
FD300K65	300K starter kit, orange	4/pkg

Jumbosep Centrifugal Device Membrane Inserts

Part Number	Description	Pkg
OD003C65	3K membrane insert, gray	12/pkg
OD010C65	10K membrane insert, blue	12/pkg
OD030C65	30K membrane insert, red	12/pkg
OD100C65	100K membrane insert, clear	12/pkg
OD300C65	300K membrane insert, orange	12/pkg

Accessories and Replacement Parts

Part Number	Description	Pkg
FD001X65	Filtrate receiver and cap	12/pkg
FD002X65	Sample reservoir and cap	12/pkg
FD003X65	Insert release	24/pkg