

## Kalrez O-Ring Kits for Sampling Accessories

Description	Part No.
Nebulizer, End Cap, and Burner Head O-Ring Kit for Plastic Nebulizer for PinAAcle AA. Includes: <b>09200272</b> (1), <b>09200275</b> (4), <b>09200276</b> (2), <b>09200277</b> (1), <b>09200278</b> (1), <b>09200369</b> (2)	<b>N3160157</b>
Nebulizer, End Cap, and Burner Head O-Ring Kit for Stainless Steel Nebulizer for PinAAcle AA. Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1), <b>09200275</b> (4), <b>09200276</b> (2), <b>09200277</b> (1), <b>09200278</b> (1)	<b>N3160156</b>
Nebulizer O-Ring Kit for Stainless Steel Nebulizer for PinAAcle AA. Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1), <b>09200276</b> (2)	<b>N3160192</b>
Nebulizer O-Ring Kit, fits Stainless Steel Nebulizers <b>03030358</b> and <b>03030352</b> . Includes: <b>09200271</b> (1), <b>09200273</b> (2), <b>09200274</b> (1)	<b>N2010247</b>
Nebulizer O-Ring Kit for Plastic Nebulizer for PinAAcle AA. Includes: <b>09200272</b> (1), <b>09200276</b> (2), <b>09200369</b> (2)	<b>N3160193</b>
Nebulizer O-Ring Kit for AA200/400 Nebulizer <b>N2010250</b> . Includes: <b>09200272</b> (1), <b>09200276</b> (2), <b>09200369</b> (2)	<b>N2010246</b>
End Cap O-Ring Kit for PinAAcle AA. Includes: <b>09200275</b> (4), <b>09200277</b> (1)	<b>N3160194</b>
End Cap O-Ring Kit for AAnalyst 700/800 End Cap, <b>00570984</b> . Includes: <b>09200277</b> (1), <b>09200341</b> (2), <b>09200342</b> (2)	<b>N2010245</b>

## PinAAcle Nebulizers and Assemblies Quick Reference

Description	Type	Needle	O-Ring Type	Complete Assembly	End Cap Assembly	Replacement Needle Assembly
Plastic High Sensitivity (standard)	Aqueous	Platinum/Iridium	Viton	<b>N3160144</b>	<b>N3160511</b>	<b>B3150525</b>
Plastic High Sensitivity	Aqueous or Organics	Platinum/Iridium	Kalrez	<b>N3160159</b>	<b>N3160189</b>	<b>N3160154</b>
Plastic High Sensitivity (suitable for aqua regia in addition to organics)	Aqueous or Organics	Tantalum	Kalrez	<b>N3160152</b>	<b>N3160189</b>	
Stainless Steel (standard)	Aqueous	Stainless	Viton	<b>N3160143</b>	<b>N3160505</b>	<b>03030354</b>
Stainless Steel	Aqueous or Organics	Stainless	Kalrez	<b>N3160158</b>	<b>N3160188</b>	<b>N3160153</b>
Stainless Steel	Aqueous	Plastic	Viton		<b>N3160505</b>	<b>03030402</b>
Stainless Steel	Aqueous	Platinum Alloy	Viton		<b>N3160505</b>	<b>03030133</b>
Description	Type			Burner Drain Assembly	Drain Bottle	
Drain Assembly	Aqueous			<b>N0400058</b>	<b>09200486</b>	
Drain Assembly	Organics			<b>N3160683</b>	<b>09200486</b>	

## AAAnalyst Nebulizers and Assemblies Quick Reference

Description	Type	Needle	O-Ring Type	Complete Assembly	End Cap Assembly	Replacement Needle Assembly
High-Sensitivity Plastic Body	Aqueous or Organics	Tantalum	Kalrez	<b>N2010250</b>		<b>N3160146</b>
Plastic Body	Aqueous	Platinum/Iridium	Viton	<b>N3150188</b>	<b>N0370392</b>	<b>N3150525</b>
High-Sensitivity Plastic Body	Aqueous or Organic	Platinum/Iridium	Kalrez	<b>B0505480</b>	<b>N0370392</b>	<b>B0501022</b>
High-Sensitivity Plastic Body	Aqueous or Organic	Tantalum	Kalrez	<b>B0505590</b>	<b>N0370392</b>	
Stainless Steel Body	Aqueous	Stainless Steel	Viton	<b>03030352</b>	<b>00570984</b>	<b>03030354</b>
Corrosion Resistant Metal Body	Aqueous or Organic	Tantalum	Kalrez	<b>03030404</b>	<b>00570984</b>	<b>03030402</b>
Platinum Alloy Nebulizer	Aqueous	Platinum/Iridium	Kalrez	<b>03030299</b>	<b>00570984</b>	<b>03030433</b>
Universal GemTip	Aqueous	Plastic w/Sapphire Tip	Viton	<b>N0370394</b>	<b>N0370392</b>	<b>N0371520</b>
High-Sensitivity GemTip	Aqueous or Organic	Plastic w/Sapphire Tip	Kalrez	<b>N0370393</b>	<b>N0370392</b>	<b>N0371522</b>

## Nebulizer Accessories

Description	Part No.
Capillary Tubing	<b>09908265</b>
Nebulizer Cleaning Wires Pkg. 5	<b>03030135</b>
Nebulizer Grease	<b>03030405</b>



Nebulizer Capillary Tubing



Nebulizer Cleaning Wires



Nebulizer Grease