

## Miscellaneous Accessories

### Injector/Detector Adapter 1/4 in.

Adapter fits on injector and detector outlet (inside oven) for use with 1/4 in. columns.

Description	Part No.
Injector/Detector Adapter 1/4 in.	<b>00080100</b>

### Liner Removal Tool

Allows easy removal of glass liners from injector port.

Description	Part No.
Liner Removal Tool	<b>N6100102</b>

### Silanized Glass Wool

Packing for wide-bore glass liners.



Description	Part No.
Silanized Glass Wool (2 oz.)	<b>00300905</b>

## Wide-Bore Adapter Kit

### 0.53 Capillary Column Adapter Kit

The AutoSystem Wide-Bore Adapter Kit allows conversion of the packed inlet for 0.53 mm capillary column operation. The kit includes 0 – 20 mL flow control element, wide-bore adapter with 1/16 in. fitting, wide-bore glass liner, and column support hanger. Can be used in both on-column and off-column modes of operation.

Description	Part No.
0.53 Capillary Column Adapter Kit	<b>N6120001</b>

## Convenient GC/MS Consumable Kits\*

Description	Part No.		
GC/MS PSS Injector Starter Kit	<b>N6100447</b>		
Contents	Pkg.	Qty.	Part No.
5.0 µL Autosampler Syringe		1	<b>N6101390</b>
Vial Locator (dongle)		2	<b>N6101182</b>
PSS Injector Viton O-Rings (250 °C)	10	1	<b>N6101747</b>
PSS Injector Kelrez O-Rings (350 °C)	1	1	<b>09921004</b>
PSS Split/Splitless Injector, 2 mm, No Wool		1	<b>N6121004</b>
Graphite/Vespel Ferrules, for 0.25 mm Columns	10	2	<b>09920104</b>
PerkinElmer Green Speta (50 pieces)		1	<b>N6621028</b>
Marathon Filament		1	<b>N6470012</b>
Aluminum Oxide Powder (3 oz.)		1	<b>04190197</b>

Description	Part No.		
GC/MS CAP Injector Starter Kit for Clarus 680, 580, 480	<b>N6100448*</b>		
Contents	Pkg.	Qty.	Part No.
5.0 µL Autosampler Syringe		1	<b>N6101390</b>
Vial Locator (Dongle)		2	<b>N6101182</b>
CAP Injector Viton O-Rings (250 °C)		10	<b>N9302783</b>
CAP Injector Kelrez O-Rings (350 °C)		10	<b>N9302782</b>
CAP Split/Splitless Injector, 4 mm, No Wool		2	<b>N6121002</b>
Graphite/Vespel Ferrules, for 0.25 mm Columns	10	2	<b>09920104</b>
PerkinElmer Green Speta (50 pieces)		1	<b>N6621028</b>
Marathon Filament		1	<b>N6470012</b>
Aluminum Oxide Powder (3 oz.)		1	<b>04190197</b>

\* **N6100448** Not compatible with Clarus 590 and 690 models  
Refer to the Clarus 590/690 consumable reference guide for more details