



## Pre-Assembled Magnetic Crimp and Screw Caps and Septa

Septa Type	Closure	Magnetic Caps (Pre-Assembled)	Part No. Pkg. 100	Part No. Pkg. 1000
Butyl (gray)	Crimp	Steel Magnetic Cap	<b>N6356561</b>	
PTFE/Butyl (red), 3.0 mm thick	Screw	Steel Magnetic Cap	<b>N9306076</b>	
PTFE/Butyl (red), 1.5 mm thick	Crimp	Steel Magnetic Cap	<b>N6356560</b>	
PTFE/Butyl (red), 1.5 mm thick	Screw	Steel Magnetic Cap	<b>N6356477</b>	
PTFE/Butyl (blue)	Crimp	Bi-Metal Magnetic Cap	<b>N6356565</b>	
PTFE/Butyl (Pharma Fix Septa)	Crimp	Steel Magnetic Cap	<b>N6356562</b>	
PTFE/Silicone (red)	Screw	Steel Magnetic Cap	<b>N6356474</b>	
PTFE/Silicone (blue/white) 0.060 in. thick	Screw	Steel Magnetic Cap	<b>N6356476</b>	
PTFE/Silicone (white)	Screw	Steel Magnetic Cap	<b>N9306077</b>	
PTFE/Silicone (natural)	Crimp	Steel Magnetic Cap	<b>N6356558</b>	
PTFE/Silicone (blue)	Crimp	Steel Magnetic Cap	<b>N6356559</b>	
PTFE/Silicone (blue)	Screw	Steel Magnetic Cap	<b>N6356475</b>	
PTFE/Silicone (blue)	Crimp	SPME Liner Steel Magnetic Cap	<b>N6356564</b>	
PTFE/Silicone (blue)	Crimp	Bi-Metal Magnetic Cap	<b>N6356566</b>	
Aluminum/Silicone	Crimp	Steel Magnetic Cap	<b>N6356563</b>	
<b>NEW</b> Ultra Low Bleed PTFE/Silicone (blue/white) 1.5 mm thick	Crimp	Bi-Metal Cap	<b>N9304181</b>	<b>N9304182</b>
<b>NEW</b> Ultra Low Bleed PTFE/Silicone (blue/white) 1.5 mm thick	Screw	Aluminium Magnetic Cap	<b>N9304175</b>	<b>N9304176</b>
<b>NEW</b> Ultra Low Bleed PTFE/Silicone (red/white) 1.3 mm thick	Screw	Aluminium Magnetic Cap	<b>N9304177</b>	<b>N9304178</b>
<b>NEW</b> Ultra Low Bleed PTFE/Silicone (white/blue) 1.3 mm thick	Screw	Aluminium Magnetic Cap	<b>N9304179</b>	<b>N9304180</b>

## Un-Assembled Aluminum Crimp Caps and Septa

Septa Type	Aluminum Caps (Un-Assembled)	Crimp Top	Crimp Top
		Part No. Pkg. 100	Part No. Pkg. 1000
Butyl (red)	Aluminum Cap, Star Spring and Septa	<b>B0159356</b>	<b>B0159357</b>
Butyl (gray)	Aluminum Cap, Star Spring and Septa	<b>N9306270</b>	<b>N9306271</b>
Butyl (gray) Stopper	Aluminum Cap (Ridged)		<b>N9303350</b>
Butyl (gray) Stopper	Aluminum Cap		<b>B0110728</b>
PTFE/Butyl (red)	Aluminum Cap, Star Spring and Septa	<b>B0104239</b>	<b>B0104240</b>
PTFE/Butyl (gray)	Aluminum Cap, Star Spring and Septa	<b>N9306266</b>	<b>N9306267</b>
PTFE/Red Rubber (red)	Aluminum Cap, Star Spring and Septa	<b>N9302979</b>	<b>N9302980</b>
PTFE/Silicone (white)	Aluminum Cap, Star Spring and Septa	<b>B0104241</b>	<b>B0104242</b>
Aluminum/Silicone	Aluminum Cap, Star Spring and Septa	<b>B0104243</b>	<b>B0104244</b>
	Aluminum Cap	<b>B0099814</b>	<b>N9302969</b>