


## Peristaltic Pump Tubing

PerkinElmer offers a wide range of high-quality peristaltic pump tubing for all of your AA, ICP-OES and ICP-MS applications. Remember to replace your tubing frequently and always keep spares on hand.



### Two-Stop Peristaltic Pump Tubing for Optima, Avio and NexION 300/350

ID (mm)	Stop Colors	Standard PVC (Pkg. 12)	Flared PVC (Pkg. 12)	Solvent Flex (Pkg. 12)	Viton® (Pkg. 12)	Santoprene® (Pkg. 12)	Silicone (Pkg. 6)
0.19	orange red	N0695476	N0773111				
0.25	orange blue	N0773117	N0773112				
0.38	orange green	N0777110	N0777042	N0777874			
0.44	green yellow		N0773113				
0.51	orange yellow		N077476		N0777446		
0.64	orange white		N0777711				
0.76	black black	09908587	N0777043	00473550	N0773118		00473552
1.02	white white				N0777442		N0777441
1.14	red red	09908585		09923037	N0773115	N0773119	N0691595
1.3	gray gray				N0777445	N0777444	N0777443
1.42	yellow yellow	N0777521					
1.52	yellow blue	09923536					
1.65	blue blue				N0777447		
1.85	green green	N0777518				N0773116	
2.06	purple purple	N0777561					
2.79	purple white			02508080			09923448
3.18	black white	N8122012					

### Three-Stop Peristaltic Pump Tubing for FIAS, FIMS and S10

ID (mm)	Stop Colors	Standard PVC Clear (Pkg. 12)	Pharmed Tan (Pkg. 12)	Solvent Resistant Yellow (Pkg. 12)
1.65	black black black	B0506058	N0773111	
1.75	white white white			B0507692
1.85	red red red	B0193160	*B3140730	
2.06	yellow blue yellow	B0193161		
2.2	purple purple purple	B0199034		
2.54	purple white purple		*B3140721	
2.79	black white black	B0508310		
0.51	orange yellow orange	B0506737		

\* Used for S10 and AS9X autosamplers



### Peristaltic Pump Tubing Starter Kit for NexION 1000/2000 – for Standard Pump

Peripump Tubing Kit (Starter Pack of Tubing Includes 4 of each commonly used sizes)				Part No. N8150060
ID (mm)	Stop Colors	Qty.	Type	
0.19	orange red	4	PVC Pure	
0.44	green yellow	4	PVC Pure	
0.51	yellow orange	4	PVC Pure	
1.3	gray gray	4	Santoprene	

### Peristaltic Pump Tubing Starter Kit for NexION 2000 P Configuration – for MP<sup>2</sup> Pump

Peripump Tubing Kit (Starter Pack of Tubing, Includes 4 of each commonly used sizes)				Part No. N8152380
ID (mm)	Stop Colors	Qty.	Type	
1.30	gray gray	4	Santoprene	
0.76	black black	4	PVC Flared	
0.51	orange yellow	4	PVC Flared	
0.38	orange green	4	PVC Flared	
0.19	orange red	4	PVC Flared	