

Flow Injection Sampling Supplies

Pump Magazines

Description	Part No.
Magazine for Peristaltic Pump One wide (>2 mm) tubing channel	B0507792
Magazine for Peristaltic Pump Two narrow (<2 mm) tubing channels	B0507791
Pump Roller Assembly	B3000251

Alternate Reductant Kit

Permits use of the FIAS-100, FIAS-400, and FIAS-400MS with an alternate reductant (e.g., SnCl₂ instead of NaBH₄).

Description	Part No.
Alternate Reductant Kit	B0508594
Includes:	
Mixing/Separation Assembly	B0507957
700 mm Single Flanged PTFE Tubing	B0191059
PTFE Tubing (110 mm)	B0191058
Red/Red Pump Tubing (Pkg. 12)	B0193160
Adapter A	B0193342
PTFE Tubing (300 mm)	B0198097

FIAS Accessories Case

Includes the most frequently used FIAS supplies, including adapters, pump tubings, connectors, sample loops, and PTFE tubings.

Description	Part No.
FIAS Accessories Case	B3000240

Flangeless Ferrule

For 1/16 in. tubing, Blue Tefzel®.

Description	Part No.
Flangeless Ferrule	09903799

Nafion® Single Tube Dryer

Continuously and highly selectively removes water from gaseous samples without affecting analytes. (Requires separate flow controller.)

Description	Part No.
Nafion® Single Tube Dryer	N9306016



Switching Valves

Description	Part No.
FI Switching Valve*, 5-Port For on-line pre-concentration. Included with all FIAS systems	B0507914
FI Switching Valve*, 8-Port For FIAS-MS	B0509095
Push-Button Assembly For manual program actuation with the FIAS-100 and FIAS-400. Included with the FIAS-100	B0505650
Disassembly Tool Valve dismantling tool	B0509590
One-Way Valve	B0198111

* For FI valve repair or replacement parts, please contact your local PerkinElmer Service office.

Electrically Heated Mantle (Cell) Adapter Cables

Description	Part No.
Cable for Adapting New RoHS FIAS/FIMS System to a Non-RoHS Heating Mantle	N2010255
Cable for Adapting New RoHS Heating Mantle to a Non-RoHS FIAS/FIMS System	N2010262