


## Nebulizer Assemblies and Spares

### Metal Body Nebulizer Parts for PinAAcle AA Instruments

For Nebulizers with Part Nos. **N3160143** and **N3160158**

Description	Part No.
<b>1</b> O-Ring sealing Nebulizer to End Cap	<b>09902239</b>
Solvent Resistant O-Ring Sealing Nebulizer to End Cap	<b>09200274</b>
<b>2</b> Cap	<b>N3161057</b>
<b>3</b> Venturi	<b>03031810</b>
<b>4</b> O-Ring for Venturi	<b>09902005</b>
Solvent Resistant O-Ring for Venturi	<b>09200273</b>
<b>5</b> O-Rings for Nebulizer Body (Each nebulizer requires a quantity of 2)	<b>09902015</b>
Solvent Resistant O-Rings for Nebulizer Body (Each nebulizer requires a quantity of 2)	<b>09200276</b>
<b>6</b> Nebulizer Body	<b>N3160133</b>
<b>7</b> Nebulizer Insert	<b>03031812</b>
<b>8</b> O-Ring for Nebulizer Insert	<b>09902047</b>
Solvent Resistant O-Ring for Nebulizer Insert	<b>09200273</b>
<b>9</b> Locking Ring	<b>N3161006</b>
<b>10</b> Spring	<b>03031806</b>
<b>11</b> O-Ring for Capillary Needle Assembly	<b>09902102</b>
Solvent Resistant O-Ring for Capillary Needle Assy.	<b>09200271</b>
<b>12</b> Capillary Needle Assembly W/O-Ring	<b>03030354</b>
Solvent Resistant Capillary Needle Assembly W/O-Ring	<b>N3160153</b>
<b>13</b> PTFE Washer	<b>03031802</b>
<b>14</b> Fine Adjustment Nut	<b>03031460</b>
<b>15</b> Sample Capillary	<b>09908265</b>




PinAAcle Stainless Steel Nebulizer  
Part No. **N3160143**

### Plastic Body Nebulizer Parts for PinAAcle AA Instrument

For Nebulizers with Part Nos. **N3160144**, **N3160152** and **N3160159**

Description	Part No.
<b>1</b> Impact Bead	<b>B0505086</b>
<b>2</b> Spacer	<b>B3150531</b>
<b>3</b> O-Rings for Nebulizer Body (Each nebulizer requires a quantity of 2)	<b>09902015</b>
Solvent Resistant O-Rings for Nebulizer Body (Each nebulizer requires a quantity of 2)	<b>09200276</b>
<b>4</b> Plastic Nebulizer Body	<b>B3161049</b>
<b>5</b> Spring	<b>09920483</b>
<b>6</b> Locking Ring	<b>N3151205</b>
<b>7</b> O-Ring for Capillary Needle Assembly	<b>09926127</b>
Solvent Resistant O-Ring for Capillary Needle Assembly	<b>09200272</b>
<b>8</b> Pt/Ir Capillary Needle Assembly W/Spring & O-Ring	<b>B3150525</b>
Solvent Resistant Pt/Ir Capillary Needle Assembly W/ Spring & O-Ring	<b>N3160154</b>
<b>9</b> Fine Adjustment Nut	<b>B3150462</b>
<b>10</b> Sample Tubing	<b>B3150458</b>



PinAAcle High-Sensitivity Nebulizer  
Part No. **N3160144**