


## Sample Introduction Kits



Part No. **N077707**



Part No. **N0791352**



Part No. **N0780546**



Part No. **N0790384**

### Cyclonic/Concentric Kits

ICP Model	Part No.
<b>Avio 200/500</b> Includes:	<b>N0790608</b>
1 Slot Quartz Torch	<b>N0790131</b>
2.0 mm ID Alumina Injector	<b>N0791183</b>
Baffled Cyclonic Spray Chamber	<b>N0791352</b>
Meinhard® Type K1 Nebulizer	<b>N077707</b>
Injector Adapter	<b>N0790141</b>
<b>Optima 8x00</b> Includes:	<b>N0780617</b>
1 Slot Quartz Torch	<b>N0780130</b>
2.0 mm ID Alumina Injector	<b>N0781014</b>
Glass Cyclonic Spray Chamber	<b>N0776052</b>
Meinhard® Type K1 Nebulizer	<b>N077707</b>
Quartz Torch Copper Ignitor Tape – Qty. 50 pcs.	<b>N0781097</b>
<b>Optima 2x00/4x00/5x00/7x00 DV</b> Includes:	<b>N0770730</b>
1 Slot Quartz Torch	<b>N0770338</b>
Quartz Bonnet	<b>N0775289</b>
2.0 mm ID Alumina Injector	<b>N0775177</b>
Glass Cyclonic Spray Chamber	<b>N0776052</b>
Meinhard® Type C High Dissolved Solids Nebulizer	<b>00472022</b>
Quartz Torch Copper Ignitor Tape – Qty. 120 pcs.	<b>N0775297</b>

### Scott/CrossFlow II Kits

ICP Model	Part No.
<b>Avio 200/500</b> Includes:	<b>N0790609</b>
1 Slot Quartz Torch	<b>N0790131</b>
2.0 mm ID Alumina Injector	<b>N0791183</b>
Scott Spray Chamber	<b>N0790384</b>
GemTip™ Cross-Flow II Nebulizer	<b>N0780546</b>
Injector Adapter	<b>N0790156</b>
Scott Spray Chamber Clamp Adapter	<b>N0790296</b>
<b>Optima 8x00</b> Includes:	<b>N0780616</b>
1 Slot Quartz Torch	<b>N0780130</b>
2.0 mm ID Alumina Injector	<b>N0781014</b>
Scott Spray Chamber	<b>N0770357</b>
GemTip™ Cross-Flow II Nebulizer	<b>N0780546</b>
Quartz Torch Copper Ignitor Tape – Qty. 50 pcs.	<b>N0781097</b>
<b>Optima 2x00/4x00/5x00/7x00 DV</b> Includes:	<b>N0770731</b>
1 Slot Quartz Torch	<b>N0770338</b>
Quartz Bonnet	<b>N0775289</b>
2.0 mm ID Alumina Injector	<b>N0775177</b>
Scott Spray Chamber	<b>N0770357</b>
GemTip™ Cross-Flow II Nebulizer	<b>N0770546</b>
Quartz Torch Copper Ignitor Tape – Qty. 120 pcs.	<b>N0775297</b>