

## FX-15, FX-20 UHPLC Pump Accessories

Description	Part No.
Replacement 0.5 µm Filter Frits for FX-15 (pkg. 4)	<b>N2911224</b>
High Pressure FX-15 Piston Seal Kit (Need 4 per FX-15)	<b>N2911221</b>
Piston Wash Bottle Kit	<b>N2601616</b>
2,500 µL Sample Syringe Starting Kit	<b>N2930313</b>

## FX-10, FX-15 and FX-20 UHPLC System Accessories

Description	Color	Size (cm)	Part No.
<b>ID PEEK Tubing</b>			
0.004 in.	Black	50	<b>N2916200</b>
0.005 in.	Red	50	<b>N2916056</b>
0.004 in.	Black	76	<b>N2916260</b>
0.005 in.	Red	66	<b>N2916262</b>
<b>Connector Tubing, Stainless Steel</b>			
Flexar ISO/Binary/Quaternary Pump to Autosampler			<b>N2916210</b>
Dual Mixer Connection (≤ 250 µL)			<b>N2916216</b>
Dual Mixer Connection (≤ 350 µL)			<b>N2916217</b>
FX-10 Pump A (Upper) to T-mixer			<b>N2916213</b>
FX-10 Pump B (Lower) to T-mixer			<b>N2916214</b>
FX-10 Purge Valve to Filter/Scavenger			<b>N2916212</b>
FX-10 Injector to Column			<b>N2916215</b>
FX-10 Pump A (Upper) to Dual Mixers			<b>N2916218</b>
FX-10 Pump B (Lower) to Dual Mixers			<b>N2916219</b>
FX-15 Pump to In-line Mixer			<b>N2916222</b>
FX-15 Pump to Dual Mixers			<b>N2916224</b>
FX-15 Injector to Column, Standard			<b>N2916223</b>
FX-15 Injector to Column, 33 cm			<b>N2916201</b>
3 cm Scavenger Cartridge Holder			<b>02580178</b>
Spheri5 Scavenger Columns (pkg. 5)			<b>02580202</b>
30 mL Priming Syringe			<b>09904849</b>
FX-15 C18 Scavenger Column Kit			<b>N2910765</b>
FX-15 C18 Scavenger Columns (pkg. 2)			<b>02580223</b>

# Maintenance Kits and Tools

Description	Part No.
Biocompatible Piston Seal Replacement Kit Includes: four seals, backup rings, and O-Rings	<b>N2910385</b>
Piston Seal Replacement Kit Includes: four seals, backup rings, and O-Rings	<b>N2910383</b>
Pulse Compensator Repair Kit Includes: diaphragm, elastomer plug, and seal	<b>N2600313</b>
Series 200 Bio LC Pump Maintenance Kit Includes: fuses, seals, O-Rings, and seal tools	<b>N2910346</b>
Series 200 Pump Maintenance Kit Includes: fuses, seals, O-Rings, and seal tools	<b>N2910345</b>
Micro pump Piston Seal Replacement Kit	<b>N2910384</b>
Check Valve Torque Wrench	<b>02540871</b>
Seal Removal Tool	<b>N2601295</b>
Insertion Tool	<b>N2601503</b>

# Manuals and Accessories

Description	Part No.
Backpressure Regulator (to generate backpressure for applications less than 250 psi)	<b>02540272</b>
Pump Priming Syringe (30 mL, plastic)	<b>09904849</b>
Series 200 Operator's Manual	<b>09938910</b>
Series 200 Pump Injector Bracket	<b>N2910349</b>