


## Plastic Body Nebulizer Parts for AAnalyst AA Instruments

AAnalyst 100/300/200/400/700/800\* High-Sensitivity Nebulizer Parts

For Nebulizers with Part Nos. **B0505480** and **B0505590**

Description	Part No.
<b>1</b> Impact Bead	<b>B0505086</b>
<b>2</b> Solvent Resistant O-Ring for Nebulizer Body	<b>B0174378</b>
<b>3</b> Nebulizer Body	<b>B0506666</b>
<b>4</b> Spring	<b>B0165332</b>
<b>5</b> Solvent Resistant O-Ring for Capillary Assembly	<b>B0190672</b>
<b>6</b> Pt/Ir Capillary Assembly	<b>B0501022</b>
<b>7</b> Locking Ring	<b>B0501025</b>
<b>8</b> Fine Adjustment Nut	<b>B0501026</b>

\*This nebulizer requires the use of end cap assembly



## AAnalyst 200/400/700\*/800\* High-Sensitivity Nebulizer Parts

For Nebulizers with Part Nos. **N3150188** and **N2010250**

Description	Part No.
<b>1</b> Impact Bead	<b>B0505086</b>
Nebulizer Spacer (not pictured)	<b>B3150531</b>
<b>2</b> O-Rings for Nebulizer Body (Each nebulizer requires a Qty. of 2)	<b>09902015</b>
Solvent Resistant O-Rings for Nebulizer Body (Each nebulizer requires a Qty. of 2)	<b>09200276</b>
<b>3</b> Nebulized Body	—
<b>4</b> Spring	<b>09920483</b>
<b>5</b> O-Ring for Capillary Assembly	<b>09926127</b>
Solvent Resistant O-Ring for Capillary Assembly	<b>09200272</b>
<b>6</b> Pt/Ir Capillary Assembly w/O-Ring	<b>B3150525</b>
Tantalum Capillary Assembly w/Solvent Resistant O-Ring and Spring	<b>N3160146</b>
<b>7</b> Locking Ring	<b>N3151205</b>
<b>8</b> Fine Adjustment Nut	<b>B3150462</b>
N/A O-Ring for Sample Tube (Each sample tube requires a Qty. 2)	<b>B0507708</b>
N/A Solvent Resistant O-Ring for Sample Tube (Each sample tube a Qty. 2)	<b>09200369</b>
N/A Solvent Resistant O-Ring Kit (Includes 1 x <b>09200272</b> , 2 x <b>09200276</b> and 2 x <b>09200369</b> )	<b>N2010246</b>

\*Requires End Cap Assembly **N0370392**

