



### SMARTsource Replacement Parts

| Description                     | Part No.        |
|---------------------------------|-----------------|
| <b>1</b> Source Lens #3         | <b>N6480149</b> |
| <b>2</b> Source Spring          | <b>N6480151</b> |
| <b>3</b> Source Lens #2         | <b>N6480148</b> |
| <b>4</b> Lens Insulator         | <b>N6480153</b> |
| <b>5</b> Source Lens #1         | <b>N6480147</b> |
| <b>6</b> EI Ion Volume          | <b>N6480144</b> |
| <b>7</b> Insulator (Ion Volume) | <b>N6480145</b> |
| <b>8</b> Repeller               | <b>N6480140</b> |
| <b>9</b> Insulator (Repeller)   | <b>N6480141</b> |
| <b>10</b> CI Ion Volume         | <b>N6480146</b> |
| <b>11</b> CI Ion Volume Disc    | <b>N6480154</b> |
| <b>12</b> Trap                  | <b>N6480142</b> |
| <b>13</b> Insulator (Trap)      | <b>N6480143</b> |
| <b>14</b> Marathon Filament     | <b>N6470012</b> |

### SMARTsource Kits

| Description   | Part No.        |
|---|-----------------|
| Complete CI Source<br>Consists of a fully assembled CI Source, ready to install | <b>N6480130</b> |
| Complete EI Source<br>Consists of a fully assembled EI Source, ready to install | <b>N6480132</b> |
| Insulator Replacement Kit (Numbers 2, 4, 7, 9, 13)                              | <b>N6480080</b> |
| EI Optics Replacement Kit (Numbers 1, 3, 5, 6, 8, 12)                           | <b>N6480081</b> |
| CI Optics Replacement Kit (Numbers 1, 3, 5, 8, 10, 12)                          | <b>N6480082</b> |