



Carbon Fiber Supported Autosampler Probes for ST Nebulizers


Self-aspiration rate @ 1 L/m Ar

Description	Part No.
0.25 mm ID, Capillary Length 80 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	N0777223
0.50 mm ID, Capillary Length 80 cm, Self-Aspiration Rate (@ 1 L/m Ar) 400 µL/min	N0777225
0.25 mm ID, Capillary Length 150 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	N0777826
0.50 mm ID, Capillary Length 150 cm, Self-Aspiration Rate (@ 1 L/m Ar) 100 µL/min	N0777827
Reset Probe Weight Kit. PTFE coated Stainless Steel Weight with Black PEEK Retainer Assembly	N0786043

PerkinElmer S10 Autosampler Probes for ST Nebulizers

Description	Part No.
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Female Fitting, 0.25 mm ID for ICP-MS. Requires but does not include N8145017	N8145139
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Female Fitting, 0.5 mm ID for ICP/ICP-MS. Requires but does not include N8145017	N8145140
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Male Screw Fitting, 0.55 mm ID for AA/ICP/ICPMS. For ICP/ICPMS requires but does not include N8145018* . For AA with a stainless steel nebulizer, the male connect should be removed so that tubing can be directly connected to the nebulizer, no additional tubing is required	N2024132
PTFE Flame Sampling Probe Assembly and Male Screw Fitting, 0.8 mm ID for AA. Requires but does not include B0172053*	N0777434
PTFE Flame Sampling Probe Assembly with Carbon Fiber Support and Male Screw Fitting, 1.0 mm ID for FIAS*	N2024131

*Keep spare ferrules on hand, **N8145237** Qty. 10



PFA MicroFlow Nebulizers with Integrated Capillaries for NexION 300/350

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA Nebulizer	20 µL/min	N8122350
PFA-50 PFA Nebulizer	50 µL/min	N8145102
PFA-100 PFA Nebulizer	100 µL/min	N8122351
PFA-200 PFA Nebulizer	200 µL/min	N8145103
PFA-400 PFA Nebulizer	400 µL/min	N8145104



PFA MicroFlow Nebulizers with Integrated Autosampler Probes for NexION 300/350

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA (Carbon Fiber Support)	20 µL/min, length 80 cm	N0777214
PFA-50 PFA (Carbon Fiber Support)	50 µL/min, length 80 cm	N0777215
PFA-100 PFA (Carbon Fiber Support)	100 µL/min, length 100 cm	N0777216
PFA-200 PFA (Carbon Fiber Support)	200 µL/min, length 100 cm	N0777217
PFA-400 PFA (Carbon Fiber Support)	400 µL/min, length 100 cm	N0777218
PFA-20 PFA (Ultem Support)	20 µL/min, length 100 cm	N8145141
PFA-50 PFA (Ultem Support)	50 µL/min, length 100 cm	N8145142
PFA-100 PFA (Ultem Support)	100 µL/min, length 100 cm	N8145143
PFA-200 PFA (Ultem Support)	200 µL/min, length 100 cm	N8145144
PFA-400 PFA (Ultem Support)	400 µL/min, length 100 cm	N8145275

PFA Nebulizers with Integrated Capillaries for NexION 1000/2000

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA Nebulizer	20 µL/min	N8152582
PFA-50 PFA Nebulizer	50 µL/min	N8152583
PFA-100 PFA Nebulizer	100 µL/min	N8152584
PFA-200 PFA Nebulizer	200 µL/min	N8152585
PFA-400 PFA Nebulizer	400 µL/min	N8152586

PFA Nebulizers with Integrated Autosampler Probes for NexION 1000/2000

Measured self-aspiration rate @ 1 L/m Ar

Description	Uptake Rate	Part No.
PFA-20 PFA (Carbon Fiber Support)	20 µL/min, length 80 cm	N8152587
PFA-50 PFA (Carbon Fiber Support)	50 µL/min, length 80 cm	N8152588
PFA-100 PFA (Carbon Fiber Support)	100 µL/min, length 100 cm	N8152589
PFA-200 PFA (Carbon Fiber Support)	200 µL/min, length 100 cm	N8152590
PFA-400 PFA (Carbon Fiber Support)	400 µL/min, length 100 cm	N8152591
PFA-20 PFA (Ultem Support)	20 µL/min, length 100 cm	N8152592
PFA-50 PFA (Ultem Support)	50 µL/min, length 100 cm	N8152593
PFA-100 PFA (Ultem Support)	100 µL/min, length 100 cm	N8152594
PFA-200 PFA (Ultem Support)	200 µL/min, length 100 cm	N8152420
PFA-400 PFA (Ultem Support)	400 µL/min, length 100 cm	N8152595