







































































Standard Pump Tubing – for NexION 1000/2000*

ID (mm)	Stop Colors	PVC Pure (Pkg. 12) Part No.	Solvent Flex (Pkg. 12) Part No.	Santoprene (Pkg. 12) Part No.
0.19	 	N8152401	N8152408	
0.25	 	N8152402		
0.38	 	N8152403	N8152409	
0.44	 	N8152404	N8152410	
0.51	 	N8152405		
0.64	 	N8152406	N8152411	
0.76	 	N8152407	N8152412	
1.14	 		N8152413	N8152414
1.3	 			N8152415



* Excluding NexION 2000 P Configuration

MP² Two-Stop Peristaltic Pump Tubing

ID (mm)	Stop Colors	Slope (L/min per RPM)	Non-Flared			Flared	
			PVC 2-Stop	Santoprene 2-Stop	Solva 2-Stop	PVC 2-Stop	Solva 2-Stop
0.13	 	0.6	N8145145			N8145194	N8145208
0.19	 	1.3	N8145146			N8145195	N8145209
0.27	 	2.7	N8145147			N8145196	N8145210
0.38	 	4.7	N8145148	N8145171		N8145197	N8145211
0.44	 	7.6	N8145149			N8145198	N8145212
0.51	 	9.5	N8145150			N8145199	N8145213
0.57	 	11	N8145151			N8145200	N8145214
0.64	 	14	N8145152			N8145201	N8145215
0.76	 	19	N8145153	N8145172	N8145176	N8145202	N8145216
0.89	 	24	N8145154		N8145177	N8145203	N8145217
0.95	 	28	N8145155		N8145178	N8145204	N8145218
1.02	 	31	N8145156		N8145179	N8145205	N8145219
1.09	 	33	N8145157		N8145180	N8145206	N8145220
1.14	 	35	N8145158		N8145181	N8145207	
1.22	 	46	N8145159		N8145182		
1.3	 	47	N8145160	N8145173	N8145183		
1.42	 	50	N8145161		N8145184		
1.52	 	51	N8145162	N8145174	N8145185		
1.65	 	55	N8145163		N8145186		
1.75	 	58	N8145164		N8145187		
1.85	 	61	N8145165		N8145188		
2.06	 	64	N8145166		N8145189		
2.2	 	65	N8145167		N8145190		
2.54	 	67	N8145168		N8145191		
2.79	 	69	N8145169		N8145192		
3.17	 	70	N8145170	N8145175	N8145193		