


Ultra Clean Filter Base Plates for GC/MS

Description	Qty.	Part No.
Ultra Clean Base plate 1 Position – 1/4 in. Brass	1	N9306800
Ultra Clean Base plate 1 Position – 1/8 in. Brass	1	N9306801
Ultra Clean Base plate 1 Position – 1/4 in. Stainless Steel	1	N9306802
Ultra Clean Base plate 1 Position – 1/8 in. Stainless Steel	1	N9306803
Ultra Clean Base plate 2 Position – 1/4 in. Brass	1	N9306804
Ultra Clean Base plate 2 Position – 1/8 in. Brass	1	N9306805
Ultra Clean Base plate 2 Position – 1/4 in. Stainless Steel	1	N9306806
Ultra Clean Base plate 2 Position – 1/8 in. Stainless Steel	1	N9306807
Ultra Clean Base plate 3 Position – 1/4 in. Brass	1	N9306810
Ultra Clean Base plate 3 Position – 1/8 in. Brass	1	N9306811
Ultra Clean Base plate 3 Position – 1/4 in. Stainless Steel	1	N9306812
Ultra Clean Base plate 3 Position – 1/8 in. Stainless Steel	1	N9306813



Replacement Filter Bundles for LC/MS

Description	Qty.	Part No.
Ultra Clean High Flow Hydrocarbon Filter Bundle: Includes (2) High Flow Hydrocarbon Filters	1 bundle of 2 cartridges	N9306823
Ultra Clean High Flow Moisture Filter Bundle: Includes (2) High Flow Moisture Filters	1 bundle of 2 cartridges	N9306824



Ultra Clean Replacement Individual Cartridge Filters for GC and GC/MS

Description	Qty.	Part No.
Ultra Clean Moisture Filter	1	N9306814
Ultra Clean Oxygen Filter	1	N9306815
Ultra Clean Hydrocarbon Filter	1	N9306816
Ultra Clean Combi (Hydrocarbon/Moisture) Filter	1	N9306818
Ultra Clean Triple (Oxygen/Moisture/Hydrocarbon) Filter	1	N9306819
Ultra Clean Triple (Oxygen/Moisture/Hydrocarbon) Helium-specific Filter	1	N9306820
Ultra Clean Triple (Oxygen/Moisture/Hydrocarbon) Hydrogen-specific Filter	1	N9306822

