


## Environmental


Analysis of Volatile Organic Compounds (VOCs) in Air Using US EPA Method TO-17.



Recommended Column: Elite-624, 60 m x 0.25 mm x 1.4 µm, Part No. **N9316006**

## Industrial


Fatty Acid Methyl Ester Contamination of Aviation Fuel by GC/MS.



Recommended Column: Elite-23, 50 m x 0.25 mm x 0.25 µm, Part No. **N9316506**

## Industrial


EPA Method 8015C for Diesel Range Organics.



Recommended Column: Elite-1, 15 m x 0.25 mm x 1.0 µm, Part No. **N9316006**

## Industrial

Determination of sulfur compounds in air by online TD-GC/FPD.



Recommended Column: Elite-5, 60 m x 0.32 mm x 0.25 µm, Part No. **N9316090**