


Clinical

Blood Alcohol Analysis.




Elite-BAC 1 Advantage, 10 m x 0.18 mm x 1.0 μ m, Part No. **N9315075**

- Baseline resolution
- Symmetrical peaks
- Run time under 2 min

Clinical


Blood Alcohol Analysis.



Elite-BAC 2 Advantage, 10 m x 0.18 mm x 1.0 µm, Part No. **N9315076**


Pharma

Residual solvents in pharmaceuticals by USP 467.

Elite-624, 30 m x 0.32 mm x 1.8 µm, Part No. **N9316203**

Petrochemical

Analysis of Permanent Gases and C1+C2 Hydrocarbons with the Elite-Carbon Column.



Column: Elite-Carbon, 2 m x 1.0 mm, 100/120 mesh, Part No: **N9303926**

Column requires an installation kit: for valve applications, Part No: **N9303450**;
for direct injections Part No: **N9303451**