

O-rings and sealing rings

| Description | Usage | Pack size | Part no. |
|--|--|-----------|----------|
| For Thermo Scientific TRACE 1300 Series GC | | | |
| Viton o-ring | Temperatures up to 300°C. Suitable for liners with OD of 6.3 mm | 10 | 0726532 |
| Graphite sealing ring | Temperatures up to 450°C. Suitable for all inlet liners above except 092004 and 09200401 | 10 | 0726005 |
| Graphite sealing ring | Temperatures up to 450°C. Suitable for use with liners 092004 and 09200401 | 10 | 0726006 |

SilTite® FingerTite ferrules



| Description | Column ID | Ferrule ID | Pack size | Part no. |
|--|-------------|------------|-----------|----------|
| For Thermo Scientific TRACE 1300 Series GC | | | | |
| SilTite FingerTite INJ / FID starter kit | 0.1-0.25 mm | 0.4 mm | * | 073610 |
| SilTite FingerTite capillary / FID starter kit | 0.1-0.25 mm | 0.4 mm | * | 073611 |
| SilTite FingerTite INJ / MS starter kit | 0.1-0.25 mm | 0.4 mm | * | 073612 |
| Replacement parts | | | | |
| SilTite FingerTite ferrule 0.4 mm | 0.1-0.25 mm | 0.4 mm | 10 | 073630 |
| SilTite FingerTite ferrule 0.5 mm | 0.32 mm | 0.5 mm | 10 | 073631 |
| SilTite FingerTite ferrule 0.7 mm | 0.53 mm | 0.7 mm | 10 | 073632 |
| SilTite FingerTite ferrule blanking | – | – | 2 | 073633 |
| SilTite FingerTite female nut | – | – | 5 | 073636 |
| SilTite FingerTite INJ base seal | 0.1-0.25 mm | – | 2 | 073640 |
| SilTite FingerTite MS adapter | – | – | 1 | 0736102 |
| SilTite FingerTite injector | 0.1-0.25 mm | – | 1 | 0736104 |

* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite FingerTite system. In addition there are five SilTite FingerTite nuts, ten SilTite FingerTite ferrules, and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.

Ferrules

| Instrument | Column ID | Size of nut | Ferrule ID | Pack size | Part no. |
|--|--------------|-------------|------------|-----------|----------|
| SilTite® metal ferrules | | | | | |
| GC-MS interface connection (starter kit) | 0.1-0.25 mm | – | 0.4 mm | 10* | 073450 |
| | 0.32 mm | – | 0.5 mm | 10* | 073451 |
| | 0.53 mm | – | 0.8 mm | 10* | 073452 |
| For TRACE 1300 split / splitless injectors (starter kit) | 0.1-0.25 mm | – | 0.4 mm | 10 | 073270 |
| | 0.32 mm | – | 0.5 mm | 10 | 073271 |
| | 0.45-0.53 mm | – | 0.8 mm | 10 | 073272 |
| | 1/32" | – | 0.81 mm | 10 | 073273 |
| Replacement SilTite® nuts | | | | | |
| SilTite metal nuts | – | – | – | 5 | 073224 |

* Includes ten ferrules, two SilTite nuts.