

O-rings and sealing rings

Description	Usage	Pack size	Part no.
Viton o-ring	Temperatures up to 300°C. Suitable for liners with OD of 6.3 mm	10	0726532
Graphite sealing ring	Temperatures up to 450°C. Suitable for all inlet liners above except 092004 and 09200401	10	0726005
Graphite sealing ring	Temperatures up to 450°C. Suitable for use with liners 092004 and 09200401	10	0726006

SilTite® FingerTite ferrules



Description	Column ID	Ferrule ID	Pack size	Part no.
SilTite FingerTite Agilent INJ / FID kit	0.1-0.25 mm	0.4 mm	*	073610
SilTite FingerTite Agilent capillary / FID kit	0.1-0.25 mm	0.4 mm	*	073611
SilTite FingerTite Agilent INJ / MS kit	0.1-0.25 mm	0.4 mm	*	073612
SilTite FingerTite Agilent INJ / FID kit	0.53 mm	0.7 mm	*	07361053
SilTite FingerTite Agilent injector kit	0.53 mm	0.7 mm	*	07361054
SilTite FingerTite Agilent capillary / FID kit	0.53 mm	0.7 mm	*	07361153
SilTite FingerTite Agilent capillary injector kit	0.53 mm	0.7 mm	*	07361154
Replacement parts				
SilTite FingerTite ferrule 0.4 mm	0.1-0.25 mm	0.4 mm	10	073630
SilTite FingerTite ferrule 0.5 mm	0.32 mm	0.5 mm	10	073631
SilTite FingerTite ferrule 0.7 mm	0.53 mm	0.7 mm	10	073632
SilTite FingerTite blanking ferrule	–	–	2	073633
SilTite FingerTite female nut	–	–	5	073636
SilTite FingerTite INJ base seal	0.1-0.25 mm	–	2	073640
SilTite FingerTite INJ base seal	0.53 mm	–	2	07364053
SilTite FingerTite capillary adapter	–	–	1	0736101
SilTite FingerTite MS adapter	–	–	1	0736102
SilTite FingerTite FID detector	–	–	1	0736103
SilTite FingerTite injector	0.1-0.25 mm	–	1	0736104

* Each starter kit includes all the parts necessary to convert one GC system (one injector and one detector) to the SilTite FingerTite system. In addition there are five SilTite FingerTite nuts, ten SilTite FingerTite ferrules, and a ferrule install tool which allows you to seat the ferrule in the correct position on the capillary column.