

1 mL to 100 mL PTFE tipped plunger syringes

Specifications			
Accuracy and reproducibility	±1% (dispensed volume)	Scale length	1 mL to 25 mL = 60 mm 50 mL = 84.2 mm 100 mL = 104 mm
Borosilicate glass barrel outer diameter (OD)	1 mL = 9 mm 2.5 mL = 11 mm 5 mL = 14 mm 10 mL = 18 mm 25 mL = 27 mm 50 mL = 32.8 mm 100 mL = 40.8 mm	Thread in plunger stem	6-32 UNC
		International standards traceability	✓



1 mL to 100 mL PTFE tipped plunger syringes

Part number	Part description and detail	Replacement needle	Replacement plunger
008020	1 mL fixed Luer Tip syringe with GT plunger	-	031842
008025	1 mL fixed Luer Lock syringe with GT plunger	-	0318441
008100	1 mL removable needle syringe with GT plunger and 5 cm 0.63 mm OD bevel tipped needle	039110	031842
008102	1 mL fixed needle syringe with GT plunger and 5 cm 0.72 mm OD bevel tipped needle	-	0318441
008110	1 mL syringe with removable needle valve and GT plunger	-	031842
008160	1 mL syringe with removable Luer Lock valve and GT plunger	-	031842
008420	2.5 mL fixed Luer Tip syringe with GT plunger	-	031852
008425	2.5 mL fixed Luer Lock syringe with GT plunger	-	031852
008500	2.5 mL removable needle syringe with GT plunger and 5 cm 0.63 mm OD bevel tipped needle	039110	031852
008502	2.5 mL fixed needle syringe with GT plunger and 5 cm 0.72 mm OD bevel tipped needle	-	031852
008510	2.5 mL syringe with removable needle valve and GT plunger	-	031852
008560	2.5 mL syringe with removable Luer Lock valve and GT plunger	-	031852
008700	5 mL removable needle syringe with GT plunger and 5 cm 0.63 mm OD bevel tipped needle	031516	031856
008760	5 mL removable Luer Lock syringe with GT plunger	-	031856
008762	5 mL fixed Luer Lock syringe with GT plunger	-	0318562
008770	5 mL syringe with removable Luer Lock valve and GT plunger	-	031856
008900	10 mL removable needle syringe with GT plunger and 5 cm 0.63 mm OD bevel tipped needle	031516	031862
008960	10 mL removable Luer Lock syringe with GT plunger	-	031862
008962	10 mL fixed Luer Lock syringe with GT plunger	-	031864
008970	10 mL syringe with removable Luer Lock valve and GT plunger	-	031862
009462	25 mL removable Luer Lock syringe with GT plunger	-	031870
009463	25 mL fixed Luer Lock syringe with GT plunger	-	031874
009472	25 mL syringe with removable Luer Lock valve and GT plunger	-	031870
009660	50 mL removable Luer Lock syringe with GT plunger	-	0312170
009670	50 mL syringe with removable Luer Lock valve and GT plunger	-	0312170
009760	100 mL removable Luer Lock syringe with GT plunger	-	0312176
009770	100 mL syringe with removable Luer Lock valve and GT plunger	-	0312176