

SGE Syringes | CTC Analytics

Product specifications and part numbers

Specifications			
Accuracy and reproducibility	±1% (dispensed volume) (±2% for 0.5 µL and 1 µL syringes)	Scale length	0.5 µL to 1 µL = 27.05 mm 5 µL to 10 µL = 54.1 mm 25 µL to 5 mL = 60 mm
Borosilicate glass barrel outer diameter (OD)	10 µL to 250 µL = 6.5 mm (except part numbers 003715 and 006720 = 8 mm) 500 µL = 8 mm 1 mL = 7.6 mm 2.5 mL = 9.7 mm 5 mL = 14 mm	International standards traceability	✓

Syringes listed in the following part number tables (pages 16-17) are suitable for:

- CTC Analytics CTC PAL
- CTC Analytics PAL RTC (PAL3)
- Shimadzu AOC-6000
- Shimadzu AOC-5000
- Thermo Scientific TriPlus RSH
- Agilent PAL3

