



□ Ordering Information

| Media Type | Target Microorganisms | Typical Incubation Time and Temperature | Order No. |
|--------------------------|-----------------------|--|------------------|
| TSA (Tryptic Soy Agar) | Total Count | 48 to 72 hrs (USP) or 1 to 5 days (EP) at 30 to 35°C | 14313--47----ACN |
| SDA (Sabouraud Dextrose) | Yeasts and molds | 5 to 7 days at 20 to 25°C | 14314--47----ACN |
| R2A | Total Count | 5 to 7 days at 20 to 28°C | 14322--47----ACN |

Accessories

| Description | Order No. |
|---|---------------|
| Microsart® e.jet Transfer Pump | 166MP-4 |
| Combi.jet manifold | 16848-CJ |
| Microsart® Base 47 mm | 1ZU---0002 |
| Silicone Pressure Tubing, pressure sided, 1 meter lengths | 1ZAS--0007 |
| Silicone tubing pressure-sided, 2 m | 1ZAS--0019 |
| Silicone tubing pressure-sided, 10 m | 1ZAS--0020 |
| Minisart® SRP25 vent filters | 17575-----ACK |