



■ Sterisart® NF gamma

Disposable units for sterility testing in isolators, individually sterile, double-packaged, gamma irradiated, needles made of flame-sterilizable stainless-steel, 10 units.

□ Ordering Information



| Description | Type of Sample | Type of Sample Container | Order No. |
|---|--|---|---------------|
| Sterisart® NF gamma with long dual-needle metal spike, sterile-vented | LVPs | Closed glass bottles with septum | 16466-----GBD |
| Sterisart® NF gamma with long needle and protective plate, inclusive sterile venting needle | LVPs SVPs | Open containers (i.e. glass ampoules, glass bottles), collapsible bags | 16467-----GBD |
| Sterisart® NF gamma with Luer or Luer Lock connection, inclusive, long needle and sterile venting needle | Medical devices | Tubing systems and bags with Luer or Luer Lock connectors | 16468-----GBD |
| Sterisart® NF gamma with syringe-adaptor and long dual-needle metal spike, sterile-vented | Pre-filled syringes | Syringes with and without needles | 16469-----GBD |
| Sterisart® NF gamma with two dual-needle metal spikes of different length, one is sterile-vented | Lyophilisates, soluble powders, liquid antibiotics | Closed glass vials with septum | 16475-----GBD |
| Sterisart® NF gamma with short dual-needle metal spike, sterile-vented | SVPs | Closed glass vials with septum | 16476-----GBD |
| Sterisart® NF gamma with long needle (side opening, with solid pointed tip, non-coring), protective plate, inclusive sterile venting needle | LVPs, SVPs, eye drops | Closed plastic containers vials ampoules, plastic containers of Blow-Fill-Seal fillings | 16477-----GBD |
| Sterisart® NF gamma with female Luer Lock connector | NEW Medical devices | Containers bags with Luer Lock male connectors | 16478-----GBD |