

## GC Accessories

Inlet Liners, Glass Wool, and Specialized Hand Tools: *Inlet Liners for PerkinElmer® (Clarus and AutoSystem)**Inlet Liners for PerkinElmer® (Clarus and AutoSystem)*

Deactivated (silanized) glass inlet liners.

**Inlet Liner, Split/Splitless Type, Single Taper FocusLiner™ Design (wool packed)**

L × O.D. × I.D. .... 92 mm × 6.2 mm × 4.0 mm



Cat. No.	Qty
2879101-U	1 ea
2879105-U	5 ea

**Inlet Liner, Split/Splitless Type, Straight FocusLiner™ Design (wool packed)**


L × O.D. × I.D. .... 92 mm × 6.2 mm × 4.0 mm



Cat. No.	Qty
2879201-U	1 ea
2879205-U	5 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 92 mm × 6.2 mm × 4.0 mm



Cat. No.	Qty
2878701-U	1 ea
2878705-U	5 ea

*Inlet Liners for PerkinElmer® (AutoSystem)*

Deactivated (silanized) glass inlet liners.

**Inlet Liner, Split/Splitless Type, Straight Design (wool packed)**

L × O.D. × I.D. .... 92 mm × 6.35 mm × 4.0 mm



Cat. No.	Qty
2631001	1 ea
2631005	5 ea
2631025	25 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 92 mm × 6.35 mm × 4.0 mm



Cat. No.	Qty
2630901	1 ea
2630905	5 ea
2630925	25 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 92 mm × 6.35 mm × 2.0 mm



Cat. No.	Qty
2631105	5 ea

**Inlet Liner, Direct (SPME) Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 92 mm × 6.35 mm × 0.75 mm



Cat. No.	Qty
2631205	5 ea

**Inlet Liner, for Packed GC (wool packed)**

Deactivated (silanized) glass inlet liner for packed injection ports.

- Fits PerkinElmer (AutoSystem)

L × O.D. × I.D. .... 112 mm × 6.0 mm × 3.0 mm



Cat. No.	Qty
2631705	5 ea
2631725	25 ea

**Inlet Liner, for Packed GC (unpacked)**

Deactivated (silanized) glass inlet liner for packed injection ports.

- Fits PerkinElmer (AutoSystem)

L × O.D. × I.D. .... 112 mm × 6.0 mm × 3.0 mm



Cat. No.	Qty
2631605	5 ea
2631625	25 ea

*Inlet Liners for PerkinElmer® (PSS Injector)*

Deactivated (silanized) glass inlet liners.

**Inlet Liner, Split/Splitless Type, Straight FocusLiner™ Design (wool packed)**

L × O.D. × I.D. .... 86.2 mm × 4.0 mm × 2.0 mm



Cat. No.	Qty
2878901-U	1 ea
2878905-U	5 ea

**Inlet Liner, Splitless Type, Straight Design (unpacked)**

L × O.D. × I.D. .... 86 mm × 4.0 mm × 2.0 mm



Cat. No.	Qty
2631301	1 ea
2631305	5 ea
2631325	25 ea