

OASIS MAX FOR ACIDIC COMPOUNDS

The Oasis MAX (Mixed-mode Anion eXchange) Sorbent has a tightly controlled ion-exchange capacity of 0.25 meq/g, ensuring reproducible SPE protocols for extraction of acidic compounds and metabolites from biological fluids. There are no silanol groups to complicate the retention mode or method development. This novel, water-wettable, polymeric sorbent is stable from pH 0–14, making method development simple and fast.

DID YOU KNOW...

When compared to other sample preparation techniques, SPE offers:

- Faster sample prep
- Lower cost
- Greater recoveries
- Greater accuracy
- Powerful enrichment of analytes
- Additional selectivity and specificity

OASIS WCX FOR STRONG BASIC COMPOUNDS

The Oasis WCX (Weak Cation eXchange) SPE material was developed to provide better sample preparation for strong bases and quaternary amines. The retention mechanism is mixed mode (both ion-exchange and reversed-phase), which improves retention for all types of basic analytes, especially strong bases.

Ordering Information

Oasis MAX Sample Extraction Products (Anion Exchange)

Description	Format	Particle Size	Qty.	P/N
Oasis MAX Cartridge	1 cc/10 mg	30 µm	100/box	186004649
Oasis MAX Cartridge	1 cc/30 mg	30 µm	100/box	186000366
Oasis MAX Flangeless Cartridge	1 cc/30 mg	30 µm	100/box	186001883
Oasis MAX Cartridge	3 cc/60 mg	30 µm	100/box	186000367
Oasis MAX Cartridge	3 cc/60 mg	60 µm	100/box	186000368
Oasis MAX Flangeless Cartridge	3 cc/60 mg	30 µm	100/box	186001884
Oasis MAX Cartridge	6 cc/150 mg	30 µm	30/box	186000369
Oasis MAX Cartridge	6 cc/150 mg	60 µm	30/box	186000370
Oasis MAX Cartridge	6 cc/500 mg	60 µm	30/box	186000865
Oasis MAX Plus Short Cartridge	225 mg	60 µm	50/box	186003517
Oasis MAX Vac RC Cartridge	20 cc/30 mg	30 µm	50/box	186000372
Oasis MAX Vac RC Cartridge	20 cc/60 mg	30 µm	50/box	186000371
Oasis MAX Vac RC Cartridge	20 cc/60 mg	60 µm	50/box	186000378
Oasis MAX µElution Plate	2 mg/96-well	30 µm	1/pk	186001829
Oasis MAX Plate	10 mg/96-well	30 µm	1/pk	186000375
Oasis MAX Plate	30 mg/96-well	30 µm	1/pk	186000373
Oasis MAX Plate	60 mg/96-well	30 µm	1/pk	186001256
Oasis MAX Plate	60 mg/96-well	60 µm	1/pk	186001205

Ordering Information

Oasis WCX Sample Extraction Products (Weak Cation Exchange)

Description	Format	Particle Size	Qty.	P/N
Oasis WCX Cartridge	1 cc/10 mg	30 µm	100/box	186004650
Oasis WCX Cartridge	1 cc/30 mg	30 µm	100/box	186002494
Oasis WCX Cartridge	3 cc/60 mg	30 µm	100/box	186002495
Oasis WCX Cartridge	6 cc/150 mg	30 µm	30/box	186002498
Oasis WCX Cartridge	1 cc/30 mg	60 µm	100/box	186002496
Oasis WCX Cartridge	3 cc/60 mg	60 µm	100/box	186002497
Oasis WCX Cartridge	6 cc/500 mg	60 µm	30/box	186004646
Oasis WCX Plus Short Cartridge	225 mg	60 µm	50/box	186003518
Oasis WCX µElution Plate	2 mg/96-well	30 µm	1/pk	186002499
Oasis WCX 96-well Plate	10 mg/96-well	30 µm	1/pk	186002501
Oasis WCX 96-well Plate	30 mg/96-well	30 µm	1/pk	186002503