

XBridge Sentry Guard Cartridges

	Dimension	P/N (2/pk)	Dimension	P/N (2/pk)
	Particle Size: 3.5 µm		Particle Size: 5 µm	
BEH C₈	3.0 × 20 mm	186003078 ²	2.1 × 10 mm	186003080 ¹
	4.6 × 20 mm	186003079 ²	3.0 × 20 mm	186003081 ²
			4.6 × 20 mm	186003082 ²
BEH C₁₈	3.0 × 20 mm	186003060 ²	2.1 × 10 mm	186003062 ¹
	4.6 × 20 mm	186003061 ²	3.0 × 20 mm	186003063 ²
			4.6 × 20 mm	186003064 ²
BEH Shield RP18	3.0 × 20 mm	186003069 ²	2.1 × 10 mm	186003071 ¹
	4.6 × 20 mm	186003070 ²	3.0 × 20 mm	186003072 ²
			4.6 × 20 mm	186003073 ²

¹Requires 2.1 × 10 mm Universal Sentry Guard Holder, p/n: WAT097958.

²Requires 3.0 × 20 mm/4.6 × 20 mm Universal Sentry Guard Holder, p/n: WAT046910.

XSelect Sentry Guard Cartridges

	Dimension	P/N (2/pk)	Dimension	P/N (2/pk)
	Particle Size: 3.5 µm		Particle Size: 5 µm	
CSH C₁₈	2.1 × 10 mm	186005252 ¹	4.6 × 20 mm	186005285 ²
	3.0 × 20 mm	186005258 ²		
	4.6 × 20 mm	186005264 ²		
HSS T3	2.1 × 10 mm	186006470 ¹	4.6 × 20 mm	186004792 ²
	3.0 × 20 mm	186004782 ²		
	4.6 × 20 mm	186004787 ²		

¹Requires 2.1 × 10 mm Universal Sentry Guard Holder, p/n: WAT097958.

²Requires 3.0 × 20 mm/4.6 × 20 mm Universal Sentry Guard Holder, p/n: WAT046910.

XTerra Sentry Guard Cartridges

	Dimension	P/N (2/pk)	Dimension	P/N (2/pk)	Dimension	P/N (2/pk)	Dimension	P/N (2/pk)		
	Particle Size: 3.5 µm		Particle Size: 5 µm		Particle Size: 3.5 µm		Particle Size: 5 µm			
MSC₁₈	3.9 × 20 mm	186000644	2.1 × 20 mm	186000652 ³	MSC₈	—	3.9 × 20 mm	186000661 ²		
	4.6 × 10 mm	186001927	3.0 × 20 mm	186000656		RP18	3.9 × 20 mm	186000646	2.1 × 20 mm	186000654 ³
			3.9 × 20 mm	186000660 ²			3.0 × 20 mm	186000658	3.9 × 20 mm	186000662 ²
			4.6 × 10 mm	186001920 ⁴						


²Requires 3.0 × 20 mm/4.6 × 20 mm Universal Sentry Guard Holder, p/n: WAT046910.

³Requires Cartridge Column Holder, p/n: 186000262.

⁴Requires In-line Guard Cartridge Holder, p/n: PSS830008.

WATERS GUARD-PAK HOLDER AND INSERTS

Waters Guard-Pak™ Holder is a compact, stand-alone housing for our unique disposable Guard-Pak Inserts. Installed in-line with your HPLC system immediately before the analytical column, the Guard-Pak Holder and inserts protect analytical LC columns against the gradual accumulation of particulates and chemical contaminants originating from the sample.



Ordering Information

Guard-Pak Holder

Description	P/N
Guard-Pak Holder	WAT088141
Guard-Pak Holder Connector	WAT080046
In-line Filters, 5/pk	WAT032472

Guard-Pak Inserts

Description	Particle size	Qty.	Pore Size	P/N
Bondapak C ₁₈	10 µm	10/pk	125Å	WAT088070 ¹
Bondapak NH ₂	10 µm	10/pk	125Å	WAT026760 ¹
Bondapak Phenyl	10 µm	10/pk	125Å	WAT026745 ¹
Nova-Pak C ₈	4 µm	10/pk	60Å	WAT015220 ¹
Nova-Pak C ₁₈	4 µm	10/pk	60Å	WAT035880 ¹
Resolve C ₁₈	10 µm	10/pk	90Å	WAT085824 ¹

¹Requires Guard-Pak Holder, p/n: WAT088141.